
iFLOW

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Table of contents

Table of contents	2
Introduction	4
General safety regulations	6
<i>Regulation regarding water and sand</i>	6
<i>Regulation regarding risk of entrapment</i>	6
<i>Regulations regarding slopes, curbs and other obstacles</i>	6
<i>Specifically for attendant support</i>	7
<i>Regulation regarding Battery and Electronics</i>	8
<i>Regulations regarding loading:</i>	8
<i>Regulations regarding the transport of a wheelchair without occupants</i>	9
<i>Regulations regarding the transport of a wheelchair with occupant</i>	10
<i>Regulation regarding EMC influences</i>	11
<i>User Instructions</i>	12
<i>Instructions regarding mounting and adjustment</i>	13
<i>Regulations regarding driving</i>	13
<i>Regulations regarding cleaning and disinfection</i>	14
<i>Storage regulations</i>	14
Used wheelchairs and the environment	16
<i>Regulation regarding the re-use of iFlow</i>	16
<i>Regulation regarding disposal of product and packaging</i>	16
Specifications	17
<i>Product specifications</i>	17
<i>Overview of product standards</i>	18
<i>Technical lifespan</i>	19

Product description	20
<i>System description</i>	20
<i>iFlow System</i>	20
<i>E-brake system and freewheel lever</i>	23
<i>Attendant support</i>	23
<i>Charger with magnetic charging plug</i>	23
Usage	24
<i>Checking the delivery</i>	24
<i>Before you drive</i>	24
<i>Using the joystick controls</i>	24
<i>Adjustability of the joystick</i>	27
<i>Using the rim support</i>	28
<i>Use of the attendant support</i>	31
<i>Battery charger indication</i>	38
<i>Wiring</i>	39
<i>Using the iFlow System Box (iBox)</i>	40
<i>Transport mode of the iFlow System Box (iBox)</i>	40
Messages and error codes	42
Troubleshooting	43
Maintenance matrix	44
Service and technical support	45
Warranty	46
<i>Manufacturer Contact Information</i>	47
Afterword	47

Introduction

Dear user,

Your wheelchair is equipped with an iFlow system. This provides additional mobility support in the form of:

- rim support iFlow Boost,
- joystick controls iFlow Joy and/or
- attendant support iFlow Care.

A combination of these features is also possible. This manual is intended to help you use the iFlow functionality correctly. This may include general explanations, safety instructions, and/or maintenance tips. Sections important to your safety are marked with a safety symbol.

If you see a **"Attention"**, **"Warning"** or **"Caution"** sign, an instruction will follow. Always follow these instructions. Failure to do so can adversely affect the operation of the iFlow system and potentially result in bodily injury or damage to the attendant support, the wheelchair, or the environment.

Warning:

- Failure to follow the safety regulations, instructions, warnings, and maintenance advice may adversely affect the operation of the attendant support and may result in bodily injury or damage to the wheelchair or the environment.
- Never change the construction of the attendant support.
- Never attempt to open the housing of the electrical components. This may cause irreparable damage. It also voids your warranty.

Attention:

- This manual is a supplement to the user manual of the wheelchair on which iFlow is mounted. All safety instructions in the wheelchair user manual also apply to a wheelchair equipped with an iFlow attendant support.
- Before using the iFlow, read this user manual and the wheelchair user manual carefully.
- Please keep this user manual in a safe place for any subsequent user.
- The iFlow system has a fully programmable electronic control system. Modifications to this programming may only be performed by Life & Mobility-authorized service technicians, dealers, occupational therapists, or product advisors.
- Use only original Life & Mobility parts and accessories.
- Never cover or remove the safety instructions and stickers.

Information for people with visual impairments can be obtained from Life & Mobility. Contact details can be found on the last page.

We'd like to keep you informed about new developments, user experiences, and other interesting information. You can regularly read the latest developments on our website www.life-mobility.com.

Kind regards,
Management & employees Life & Mobility

General safety regulations

In this chapter we would like to inform you about a number of regulations that are important for safety when using the iFlow.

Regulation regarding water and sand

- The iFlow system is splashproof (IPx4). Therefore, driving through (deep) puddles should be avoided.
- Never spray the wheelchair on which the iFlow system is mounted with a pressure washer or garden hose, as this may irreparably damage the iFlow system.
- Do not allow the iFlow system to come into contact with seawater: seawater is aggressive and will damage the iFlow system.
- Do not allow the iFlow system to come into contact with sand: sand can penetrate into the rotating parts of the iFlow, causing unnecessary rapid wear.

Regulation regarding risk of entrapment

When using the wheelchair, be careful not to get your hands caught between any parts of the wheelchair or the iFlow system. This includes the wheel spokes, handrims, the frame, the brake levers, and the leg rests. You could also get caught between the wheelchair and surrounding obstacles, such as doorframes or furniture. Therefore, always exercise caution and avoid unnecessary risks.

Regulations regarding slopes, curbs and other obstacles

Warning:

- Driving up or down a curb (at an angle) can cause serious damage to the iFlow system. Therefore, always use a level driveway when driving up or down a sidewalk.

- Before releasing the parking brakes using the freewheel levers, check that this will not cause the wheelchair to move unintentionally.

Specifically for attendant support

- This can create dangerous situations, especially on slopes. Do not engage freewheel mode unless absolutely necessary. In these situations, have the occupant exit the wheelchair first and ensure there is sufficient assistance to hold the wheelchair. (See also 'Driving without attendant support').

Attention:

- Avoid slopes steeper than 7°.
- Always keep your speed under control when descending.
- Never drive down a slope in freewheel mode.
- The iFlow system consumes much more power on inclines than on a flat road. This can significantly reduce the range. Pay attention to the battery indicator to avoid unexpected stops.
- When going downhill, the system's stopping distance may be significantly longer than on a level surface.

Regulation regarding Battery and Electronics

(See also chapter 'Use').

Regulations regarding loading:

- Use only the supplied iFlow battery charger to charge the iFlow battery. Charging batteries other than the iFlow battery may damage the system.
- Avoid short circuits between the battery charger terminals.
- Never charge the iFlow with an occupant in the wheelchair.
- The operation of the battery charger is reserved for those who can properly understand the instructions and safety precautions.
- Keep the battery charger out of reach of children.
- Never open the battery charger.
- Have maintenance carried out only by certified service technicians.

Caution:

- Stop the charging process by unplugging the power cord:
 1. In case of abnormal smell, leakage, deformation of the charger, battery or other parts.
 2. If the charging process takes longer than usual.
- The battery charger may heat up during the charging process. This is normal. Therefore, never cover the battery charger while charging.
- Before each use, check that the insulation of the battery charger plug and cables is still in good condition.
- Do not use any other cables than the supplied iFlow charging cables.

- Ensure that electronic components remain protected from excess water, e.g. during cleaning operations.
- Please refer to the section 'Disposal of the product and packaging' for information on disposal of this product.

Regulations regarding the transport of a wheelchair without occupants

There are a number of additional regulations and points of attention regarding transport without occupants.

Warning:

- When transporting the wheelchair in a taxi, ALWAYS put the parking brake on.

Attention:

- If your iFlow system is equipped with an E-brake, it will NOT automatically engage the parking brake mode when freewheel is engaged!

Attention:

- If your iFlow system is equipped with an e-brake, do not mistake the freewheel levers for the parking brake. The levers should remain pointed vertically upward.

Attention:

- Always switch off the iFlow system during transport.
- Lift the wheelchair only by the fixed frame parts and not by any parts of the iFlow such as the attendant support, joystick controls or wheels.
- Make sure controls remain clear during transport. They can be damaged by impact or pinching.

Regulations regarding the transport of a wheelchair with occupant

The iFlow system has been tested for passenger transport when combined with the Life & Mobility Roxx, Canto, and Match wheelchairs. Agreements may also exist with other wheelchairs. Contact Life & Mobility for the possibilities. Other combinations have not been tested according to this standard and are therefore not suitable for passenger transport with occupants. These wheelchairs with the iFlow system are designed to be secured in a forward-facing position in a motor vehicle. The iFlow system has been successfully tested as a crash-test safe system in accordance with ISO 7176-19:2008. A wheelchair with the iFlow system can therefore be transported as a seat with occupant in a taxi bus (see also the identification sticker on the seat). The condition of the occupant and the doctor's advice are decisive for safe transport. Life & Mobility recommends transferring the occupant to a fixed seat, if possible.

Attention:

- NEVER use the power assist to drive the chair and its occupant into a taxi van via ramps. The ramp angle is too steep for this. Use the lift instead.
- Switch off iFlow when transporting the wheelchair with occupant.
- Transport with occupants is only permitted if the wheelchair is fitted with a crash test safe sticker.
- Always follow the instructions in the wheelchair manual for safe passenger transport.

Warning:

- When transporting the wheelchair in a taxi, ALWAYS put the manual parking brake on, if available.
- If your iFlow system is equipped with an E-brake, it will automatically engage the parking brake when the system is switched off.
- Don't confuse the freewheel with the parking brake. The parking brake is automatically activated when the attendant support is switched off.

Regulation regarding EMC influences

The iFlow system is EMC (Electromagnetic Compatibility) approved according to EN-12184, which refers to EMC/EMI requirements for the use of medical devices. Tests have shown that the iFlow is not sensitive to or harmful to mobile communication equipment and electronically controlled devices.

Nevertheless, be careful when using the iFlow system simultaneously with mobile or electronic equipment. Use near strong electric or magnetic fields, such as those found in industrial environments, is also not recommended. If you notice any problems, the attendant support should be turned off.

User Instructions

The iFlow system is a modular system that can be configured to meet specific user requirements. Configured as driver support for the attendant, the system imposes no additional requirements on the occupant. The requirements for the occupant specified in the wheelchair's user manual remain in effect. When configured as driver support for the occupant (joystick control or handrim support), there are additional requirements in addition to those specified in the wheelchair's user manual:

- You should not use iFlow if you are under the influence of drugs that may impair your ability to drive.
- You must be physically and mentally able to operate iFlow in the intended environment.
- Your vision must be sufficient to operate iFlow safely in any situation that may arise during its use.
- Make sure you're always visible to other road users. As an iFlow user, you're particularly vulnerable in traffic.
- Please observe applicable traffic regulations.
- Avoid quiet routes. Make sure help can be called in an emergency.

Attention:

- In the event of a malfunction or a dead battery, the attendant support will be disabled. The attendant must be able to bring the wheelchair to a safe position under these circumstances. This can be done, for example, by placing the wheelchair in freewheel mode and pushing it further under your own power, or by leaving the brake engaged and calling for assistance.
- When using the iFlow as an attendant support on slopes, the condition of the attendant also determines the maximum slope on which the wheelchair can be used safely.

Instructions regarding mounting and adjustment

Warning:

- Before using a wheelchair with iFlow, check that all parts are properly assembled and secured.

Attention:

- When adjusting the height and angle of the push handles, check that the cables do not become pinched.

Regulations regarding driving

(See also chapter 'Use').

Warning:

- Switch off the attendant support if the wheelchair is left unattended.

Regulations regarding cleaning and disinfection

All parts of the iFlow can be cleaned with a slightly damp cloth.

Caution:

- Do not use bleach or solvents. Do not use a pressure washer to clean the iFlow. Clean the housings with a non-abrasive cleaner.

Attention:

- Clean the system parts regularly.

All iFlow components can be disinfected with a commercially available disinfectant wipe if desired. Do not wipe the components dry after disinfection; allow them to air dry. If reusing the iFlow components, they should be disinfected before being given to a new user.

Storage regulations

Wheelchairs equipped with the iFlow system cannot be reduced in size. For long-term storage, the iFlow system box (iBox) transport mode must be activated; see the chapter on using the iFlow system box (iBox) for more information. There are a few additional considerations regarding iFlow storage.

Caution:

- The wheelchair must be stored in a dry place. The humidity should not exceed 80%. Furthermore, condensation should not form on the wheelchair, for example, due to temperature fluctuations.

Attention:

- Ambient storage temperature: -10 to +60°C.
- See the maintenance matrix in this user manual for information on the minimum frequency of charging the battery.

* For long-term storage with battery: 10-30°C.

Used wheelchairs and the environment

Regulation regarding the re-use of iFlow

iFlow can be reset for a new user. This should only be performed by qualified personnel. The necessary information can be found in the service manual.

Regulation regarding disposal of product and packaging

The system box (iBox) contains the system's battery. Do not dispose of the battery with normal household waste. You can return your discarded iFlow system to your dealer. The dealer will return the iFlow system to Life & Mobility. We will ensure it is disposed of correctly.

The remaining iFlow components are considered electronic waste. Please check with your local government for applicable regulations regarding the disposal of these components. For environmental reasons, we recommend that you separate the packaging and the system as much as possible so that the components can be recycled efficiently. Your dealer can tell you where to return your discarded iFlow.

- * *The maximum value that can be set may vary per country.*
- ** *In practice, the range may be lower. This is negatively affected by the weight of the occupant, weather conditions, terrain, slopes, and battery age.*
- *** *For long-term storage with battery: 10-30°C*

Specifications

Below you will find the specifications of the iFlow system. Other specifications for the wheelchair can be found in the user manual for the specific wheelchair on which the iFlow system is installed. More information is available on our website.

Product specifications

- Intended use: Indoors and outdoors on paved surfaces
- Maximum user weight 136 kg.
- Maximum safe slope: 7 degrees
- Additional weight iFlow attendant support: 20-26 kg
- Maximum technical speed*:
 - Rim support 6 km/h
 - Joystick 10 km/h
 - Attendant support 6 km/h
- Maximum obstacle height: 15 mm
- Maximum braking distance at maximum speed (programmable): 2.1 meters
- Battery type: LiFePO4 Li Ion 25.6 V
- Battery capacity: 12 Ah
- Battery weight: 4.3 kg (complete system box (iBox))
- Charging time: approx. 4 hours for full charge
- Brushless motor: 450W
- Range: 15 km**
- Charger input: 100 – 240V 50/60Hz
- Charger output: 28.8V DC/ 3A
- Ambient temperature during use: -5°C to 45°C
- Ambient temperature during charging: 0°C to +45°C
- Ambient temperature for storage and transport: -10°C to +60°C***

The iFlow system can be used safely in combination with the following wheelchairs from Life & Mobility: Roxx, Canto, Match or additional covenants.

Component	Weight (kg)
System box (iBox)	4,3
26" wheels with electric brake (incl. freewheel levers)	18,7
26" wheels with rim sensors, without electric brake	19,7
12.5" wheels, without electric brake	13,8
12.5" wheels, with electric brake	17
Joystick controls	0,4
Rim control panel	0,1
Complete control panel attendant support	2,0

Overview of product standards

The iFlow companion support, in combination with the wheelchair, meets the following standards:

- EN 12183:2022 Manually wheelchairs - Requirements and test methods
- EN 12184:2022 Electrically powered wheelchairs, scooters and associated chargers
- ISO 7176-8:2014 Requirements and test methods for static, impact and fatigue strenght
- ISO 7176-9:2009 Climate tests for electric wheelchairs
- ISO 7176-14:2022 Power and control systems for electrically wheelchairs and scooters
- ISO 7176-19:2022 Wheelchairs for use as seats in motor vehicles
- ISO 7176-21:2009 Requirements and test methods for Electromagnetic compatibility

of electrically powered wheelchairs and motorized scooters

- UN/DOT 38.3 Transport of Dangerous Goods,
- NEN-EN-IEC 62133-2:2017 Safety requirements for portable sealed cells and for batteries for use in portable applications

Technical lifespan

Life & Mobility assigns this product a technical lifespan of 5 years. This means the product must be repairable, safe, and usable during this period.

The above statement applies only when the user uses and maintains the product in the manner prescribed by Life & Mobility, as indicated by the CE standard, and using original accessories.

Product description

This chapter describes the various components of the iFlow system. We'll first show how they work together. Then, we'll discuss the individual parts one by one.

System description

iFlow is a modular system of auxiliary drives for manual wheelchairs. It comprises three distinct systems, which can be combined thanks to their modular design. The three systems are:

- iFlow Boost -> Rim support
- iFlow Joy -> Joystick control
- iFlow Care -> Attendant support.

The following image shows a wheelchair with all three systems combined. This provides a complete overview of all the options. In practice, a choice will be made based on the indication.

iFlow System

- A. iFlow System Box (iBox) (Battery and system controller)
- B. Rim support
- C. Rim support Control module
- D. Optional on/off button for rim support
- E. Wheel with drive motor
- F. Joystick for joystick control
- G. E-brake braking system + freewheel lever
- H. Attendant support
- I. Charger with magnetic charging plug

<p>iFlow System Box (iBox)/System Controller The system box (iBox) processes the input signals and controls the motor. The lithium (LiFeP) battery is also integrated into this box. Up to three input signals can be connected.</p>	
<p>Control module for rim support This module is the interface for the rim support. It displays the battery status and configures riding profiles. It also charges the system.</p>	
<p>Optional on/off button for rim support The handrim support can optionally be switched on and off with a simple button. This can be operated by the user. In that case, the Control module can be mounted on the back/side of the chair.</p>	
<p>Rim support The rim support consists of a separate rim that can rotate relative to the motor wheel. By operating this, the system receives the input needed to control the motor.</p>	
<p>Drive motor for joystick and attendant control The iFlow system is equipped with brushless, direct-drive motors. These are integrated into the spoked wheels (20-26 inches) and the small unspoked wheel (12.5 inches). The motors are internally equipped with an additional encoder to ensure highly precise control.</p>	

<p>Joystick for joystick control The joystick provides input and serves as the interface for joystick control. It allows you to select driving profiles and operate the horn. The color display also provides access to the iFlow system settings. The charging connector is located at the rear.</p>	
<p>E-brake system and freewheel lever In the case of joystick control and optionally with attendant support, the iFlow has an electric brake. This functions as an automatic parking brake. In the event of a malfunction, the brakes can be released using the two freewheel levers. These are marked with the mandatory freewheel symbol.</p>	
<p>Attendant support The companion support system consists of an intuitive handlebar with a central interface and two controllable handles. These are equipped with sensors that provide input to the system. The lever on the right side activates the system.</p>	
<p>Charger with magnetic charging plug The iFlow system is charged using a 3A charger. This is connected to the modules with a magnetic charging connector.</p>	

Usage

This chapter contains the instructions essential for the correct use of the iFlow system. Life&Mobility is confident that you will enjoy using the system if you follow these instructions carefully. If you have any questions or concerns, please contact your supplier or Life&Mobility.

Checking the delivery

Before using the iFlow system, check that the delivery matches your order. If you notice any transport damage upon delivery, please report it to your dealer.

Before you drive

Start slowly and take your time to get used to the support.

Using the joystick controls

Control panel Joystick

The iFlow system can be equipped with a joystick control. In that case, the joystick control is attached to the wheelchair armrest using either a fixed or swing-away mount. The wheelchair occupant is the intended user of the joystick control.

- A. Battery indication
- B. Driving profile indication
- C. Driving profile buttons
- D. On/off button
- E. Horn button
- F. Joystick
- G. Charging port (bottom of Joystick)

Enabling and disabling the joystick control
Press the power button (D) to turn the system on. The joystick control display will now light up. Press the button again to turn the system off when the system is on.

Perform an emergency stop using the joystick control

If an emergency occurs during assisted driving, an emergency stop can be made by switching off the system using the on/off button (D)

Attention:

- Use the emergency stop function only in emergency situations. See the chapter "Driving with Joystick Controls" for instructions on how to stop the system during normal operation.

Using the horn

Press the horn button (E) to activate the horn. The horn volume can be adjusted by your occupational therapist or dealer. The horn can also be turned off.

Driving with the Joystick Controls

If you see a driving profile (a number of 1 or higher) you can drive using the joystick controls, see image.

- *To ride*
Push the joystick (F) forward to drive forward. Push the joystick back to drive backward. Push the joystick to one side to turn in that direction. The further the joystick is pushed from the center, the faster the acceleration and the higher the final speed you can achieve.

- **To come to a standstill**
Move the joystick toward the center position to reduce speed. Release the joystick completely to come to a stop.
- **Change driving profile**
Use the drive profile buttons (C) to change drive profiles. Press the up arrow button to set a higher profile. Press the down arrow button to set a lower profile. The iFlow joystick control can contain up to five drive profiles. Each profile can be configured with a different driving behavior. Consult your occupational therapist for the most optimal drive profiles.
- **Maneuvering mode**
Activate maneuvering mode by pressing the down arrow button when driving profile 1 is selected and the wheelchair is stationary. An "M" will appear on the display. At this point, the electric brake is deactivated and the wheelchair can be moved without assistance. Deactivate maneuvering mode by pressing the up arrow button.

Messages on the joystick control

If the joystick is operated during power-up, this symbol (Fig. A) will be displayed. The joystick is not in the neutral position, and the system will not boot further. Release the joystick to continue booting the system.

Freewheel active

The message 'Freewheel active' (fig. B) is displayed when the Electric brake is in the freewheel mode, see section 'Freewheel'.

System will soon go into sleep mode (fig. C)
The iFlow system can be set to automatically shut down after a certain period of inactivity. Consult your dealer or occupational therapist if you wish to change this setting. If this function activates and the system goes into sleep mode within 5 minutes, you will see the following message. Press iFlow to cancel.

Adjustability of the joystick

The joystick control can be attached to the wheelchair by means of a fixed mounting or a parallel swing-away mounting.

Fixed mount

A joystick control mounted on a fixed mount can be adjusted to the desired position and angle. Standard hand tools are required for this. Consult your occupational therapist, dealer, or advisor.

Parallel swing-away recording

A joystick control mounted on a swing-away mount allows the user to push it aside. It can also be adjusted to the desired position and angle. Standard hand tools are required for this. Consult your occupational therapist, dealer, or advisor.

Joystick control in combination with driver support

For a system equipped with both joystick control and attendant support, the explanation below applies. The chapter "Using Attendant Support" provides further explanation of attendant support.

A dash instead of a driving profile number on your joystick screen means you cannot drive; the attendant is using the system. The system can be configured in two ways:

- Automatic control override. You can take over control of the system once the chair has come to a stop.
- Manual control transfer by an attendant. The attendant must press the "Control Transfer" button to transfer control to the joystick.

When the message "Manoeuvring mode attendant active" appears, the attendant has activated manoeuvring mode. This means the electric brakes are not engaged, allowing the chair to be maneuvered without assistance.

Using the rim support

Control panel Rim support

- A. Battery indication
- B. Driving profile indication
- C. On/off button
- D. Driving profile selection button
- E. Charging port
- F. Grab rims with force sensors

The iFlow system can be equipped with a handrim support. In that case, the handrim Control module is attached to the chair frame. The wheels are equipped with sensors that detect how hard you press on the handrims. The wheelchair occupant is the intended user of the handrim support.

Switching the rim support on and off

Press the power button (C) to turn on the system. Several LEDs will now light up. Press the button again to turn the system off when the system is on.

Riding with the rim support

- **Drive forward**
Push both handrims (F) forward to drive forward. The harder you push, the stronger the support becomes.
- **Reversing**
Push both handrims back to reverse. The harder you push, the stronger the support becomes.
- **Coming to a stop while driving**
To come to a stop, hold the handrims against the handlebars while driving.
- **To steer**
Push harder on one rim than the other rim to make a turn.
- **Change driving profile**
Change the driving profile by pressing the driving profile button (D). This will activate a higher driving profile. If the highest driving profile is already active, the lowest driving profile will be activated. The iFlow push rim support can support a maximum of three driving profiles and a minimum of one. Each profile can be configured with a different driving behavior. Consult your occupational therapist for the optimal driving profiles.

Notifications on the rim support

A pulsating red drive profile indicator LED is displayed when the wheelchair is turned on while the handrims are not in the neutral position. When these controls are released, the pulsating red light disappears and the system continues to power up. You can drive as soon as you see a white drive profile indicator LED.

Adjustability of the rim support control module

The handrim support control panel can be adjusted in angle and position. Standard hand tools are required for this. Consult your occupational therapist, dealer, or advisor.

System with rim support and attendant support

For a system equipped with both a handrail and a companion support, the explanation below applies. Further information about the companion support is provided in the chapter "Using the companion support."

If none of the driving profile indicator LEDs light up, you cannot drive with assisted driving.

The system can be set up in one of two ways:

- Automatic control override. As soon as the wheelchair has come to a stop and one of the driving profile indicator LEDs on the handrail support lights up white, you can resume assisted propulsion.
- Manual control override by an attendant. As soon as the attendant presses the "Control Override" button and one of the driving profile indicator LEDs on the handrail support lights up white, you can resume assisted handrail operation.

Using the optional on/off button

The system can optionally be equipped with an on/off button. This is attached to the wheelchair frame with a bracket.

- A. On/off button
- B. On/off indicator light ring
- C. Angle adjustment

Turn the system on by pressing the power button. A green light ring on the power button indicates the system is on. Turn the system off by pressing the power button again. If the light ring does not light up, this indicates the system is off.

Adjusting the on/off button

The optional on/off button can be adjusted in angle and position. Standard hand tools are required for this. Consult your occupational therapist, dealer, or advisor.

Use of the attendant support

Control panel Attendant support

The iFlow system can be equipped with an attendant support. In that case, the handlebar with the display is mounted on the back of the chair. The wheelchair attendant is the intended user of the attendant support.

- A. Battery indication
- B. Driving profile indication
- C. Driving profile buttons
- D. On/off button
- E. Horn button
- F. Control surrender button
- G. Maneuver button
- H. Charging input
- I. Hand detection lever. The following icon (J) is displayed on the lever.

If equipped with a tilt lever, its operation is indicated by the following icon (K)

Enabling and disabling the attendant support
Press the power button (D) to turn on the system. The screen will now light up. Press the button again to turn off the system when the system is on.

Behavior of attendant support when iFlow system is off (system **with Electric brake)**

An iFlow system with an electric brake will bring the wheelchair to a stop if it is switched off while driving. If an emergency occurs during assisted driving, an emergency stop can be initiated by deactivating the system using the on/off button (D).

Attention:

- Use the emergency stop function only in emergency situations. See the "Driving with Attendant support" section for instructions on how to stop the system during normal use.

Behavior of attendant support when iFlow system is off (system **without Electric brake)**

An iFlow system without an electric brake functions like a manually operated wheelchair when the system is turned off. See your wheelchair's user manual for further explanation. See the "Driving with Attendant support" section of this user manual for instructions on how to stop the system when the iFlow system is turned on.

Caution:

- Use the tire brake or foot brake to park an iFlow system without an electric brake. This prevents unintentional driving or rolling away.

Horn

Press the horn button (E) to activate the horn. The horn volume can be adjusted by your occupational therapist or dealer. The horn can also be turned off.


Driving with the attendant support

- Forward/backward driving
Squeeze the detection lever (I) and push the handles in the direction you want to drive.
- To come to a standstill
To slow to a stop, stop applying force to the handles. Release the detection lever to stop quickly.
- To steer
Push the right handle harder than the left handle to make a left turn. Push the left handle harder than the right handle to make a right turn.
- Change driving profile
Use the driving profile buttons (C) to change driving profiles. Press the up arrow button to set a higher profile. Press the down arrow button to set a lower profile. iFlow's attendant support can contain up to 5 driving profiles and at least 1 driving profile. Different driving behaviors can be set for each profile. Consult your occupational therapist for the most optimal driving profiles.
- Maneuvering mode
Activate maneuver mode by pressing the Maneuver button (G) when the wheelchair is stationary. An "M" will appear on the display. At this point, the electric brake is deactivated and the wheelchair can be moved without assistance. Deactivate maneuver mode by pressing the Maneuver button (G) again.


Release push support

This symbol appears when the wheelchair is turned on while the push handles or the detection lever (I) are not in the neutral position. When these controls are released, the symbol disappears and the system continues to start.


'Freewheel active'

The message 'Freewheel active' is displayed when the Electric brake is in freewheel mode, see section 'Freewheel'.

'System will soon go into sleep mode'

The iFlow system can be set to automatically shut down after a certain period of inactivity. Consult your dealer or occupational therapist if you wish to change this setting. If this function activates and the system goes into sleep mode within 5 minutes, you will see the following message. Press iFlow to cancel.

Adjustability of the attendant support

The attendant support is mounted on a height-adjustable mounting bracket. Adjustment methods vary depending on the chair type. Consult your wheelchair's manual for the correct adjustment instructions.

System with attendant support and joystick control or rim support

For a system equipped with both a handrim support and a joystick control, the explanation below applies in addition to attendant support. Further explanation of these components is provided in the chapters on using the handrim support and joystick control.

If you see a driving profile (a number of 1 or higher) you can drive with support.

If you see a dash instead of a driving profile number, this means you cannot drive with assisted driving; the system is using the joystick control or rim support.

The system can be set up in two ways:

- Automatic steering override. You can take control of the system and drive with assistance once the seat has come to a stop.
- Manual control override by attendant. Press the "Control Override" button when the wheelchair is stationary to take control of the system and drive with assistance.
- When the message "Occupant is in control" appears, the occupant has activated maneuvering mode. This means the electric brakes are not engaged, allowing the seat to be maneuvered without assistance. Maneuvering mode does not affect your ability to drive with assistance. See the chapter "Driving with attendant supportance" for instructions on how to use maneuvering mode as an attendant.

Using the Electric Brake


If the system is equipped with joystick control, each wheel has an e-brake. If the system is equipped with attendant support and/or handrim support, the electric brake is optional.

During assisted driving, the electric brakes are electronically controlled by the iFlow system. In the event of a system malfunction, the motors can be placed in freewheel mode. To do this, operate the freewheel levers on the left and right motors. These are identified by the symbol shown here.

Read battery status

The remaining battery capacity is indicated by five blocks on the control modules. Below 20%, the leftmost block turns red. Below 10% battery capacity, the leftmost block flashes.


Charging the iFlow system

The iFlow battery charger supplied with the iFlow system is exclusively intended for charging the iFlow system.

General specifications of the iFlow battery charger:

Only applicable for:	iFlow Lithium Iron Phosphate, (LifePO) 12Ah
Temperature range:	+10-40°C, IPX1
Supply voltage:	100-240VAC
Loading capacity:	3A
Output voltage DC:	28,8 V (24 V)

If the battery needs to be charged, the included iFlow battery charger must be connected to a charging port. The charging port is located on the Control modules (joystick, rim support control module, or attendant support). The available control modules on the wheelchair depend on the selected iFlow configuration.

Each charging plug features a magnetic contact that simplifies connection. This magnetic connection ensures frequent connection and disconnection. It also reduces the risk of damage if the connection is accidentally disconnected.

The battery doesn't need to be empty before recharging. You can charge it at any time. You can also interrupt the charging process at any time. Charging a fully discharged battery takes approximately 4 hours.

To charge, do the following:

1. Connect the magnetic charging connector to a power outlet on the wheelchair.
2. Plug the charger's power cord into a wall outlet. While charging, the light on the charger will glow red. The charging animation will be displayed on the joystick control. The number of stationary blocks indicates charging progress, with each stationary block representing 20% of the stored capacity.

You cannot drive while charging. To resume

driving, do the following:

1. Unplug the appliance.
2. Disconnect the magnetic charging connector from the wheelchair.
3. Turn the system off and then on again.

Attention:

- Whether the system is on or off does not affect the charging process.
- You cannot use the system while charging.
- Charge the battery only at an ambient temperature of +10 to +40°C*.
- At ambient temperatures (during use and storage) lower than -10°C or higher than 60°C, the battery will have a shorter life and will need to be charged more often.

Battery charger indication

Lamp	Situation
Green-red alternating	Charger not connected, battery not charging
Red	Battery is charging
Green	Battery is fully charged

Auto power off

The iFlow system can be set to automatically shut off after a certain time. Optionally, the system will emit a beep before shutting itself off. This beep is disabled by default. Consult your dealer if you wish to change the setting for this function.

Motor cables

Module cables

Wiring

Motor cables

The iFlow wheels with integrated motors are connected to the iFlow System Box (iBox) via cables with connectors. Connecting and disconnecting these cables affects the functioning of the iFlow system and should therefore only be performed by a dealer or authorized service technician.

Module cables


A joystick control, attendant support, rim support Control module, and on/off button are connected to the iFlow System Box (iBox) via a module extension cable. Connecting and disconnecting these cables affects the functioning of the iFlow system and should therefore only be performed by a dealer or authorized service technician.

Using the iFlow System Box (iBox)

Each iFlow system is equipped with a system box (iBox). The system box contains the battery and the system's Control module. The iFlow System Box (iBox) is mounted to the wheelchair and cannot be removed during use.

- A. System box (iBox)
- B. Motor connections (2x)
- C. Function button for transport mode
- D. Module connections (3x)
- E. iFlow Registration/serial number sticker

Transport mode of the iFlow System Box (iBox)

The iFlow system box (iBox) can be set to transport mode if desired. The following scenarios are examples where this mode is desirable:

- long-term storage,
- transportation in an aircraft or other mode of transportation,
- service and maintenance
- error situations
- in case of fire or other emergency situations.

Enabling storage mode disconnects the battery in the iFlow System Box (iBox) from the rest of the electronics.

Transport mode is activated by pressing the function button (C) on the system box (iBox) for

System box (iBox) with indication Serial number sticker

System box (iBox) on the Canto Nxt

System box (iBox) on the Roxx

System box (iBox) on the Match

8 seconds.

Attention:

- When storage mode is enabled, the system cannot be switched on and therefore iFlow cannot be driven.
- To put a System Box (iBox) in Transport mode back into use, it must be activated. The iFlow System Box (iBox) is reactivated by connecting an iFlow charger to one of the charging inputs on the controls.

User Manual

Manufacturer → Logistiekweg 7
7007CJ Doetinchem
T +31 (0)314-328 000
@ info@life-mobility.com

Medical tool → **MD**

Product name → **REF iFlow** max. user weight 136 kg 2024-09

Serial number → **SN** (01) 08719327267301 (21) S_IF000014

CE

Date of manufacture

Type Identification Sticker

Messages and error codes

The control modules display a message in certain situations. The chapters on using the control modules cover messages that may appear during normal use. The table below shows other messages and error codes.

Indication on joystick control or attendant support	Indication of rim support	User action
System error: system box (iBox) System error: battery System error: brake left System error: brake right System error: rim support System Error: Joystick System error: attendant support System Error: Charger		Consult your dealer.
System error: left engine System error: right engine System error: left rim sensor System error: right rim sensor		Consult your dealer.
Connect 'component' (for example: Connect left brake)		Check that all cables are connected correctly or consult your dealer.
System temperature too high, let it cool down System temperature too low, let it warm up Engine limitation due to temperature		Allow the system to warm up or cool down, preferably in an environment that is at room temperature.
System not configured Device not configured System configuration conflict		Consult your dealer
Battery empty Motor limitation by voltage		Charge battery
	Blinking	
Maintenance required (Message only during start-up, disappears as soon as the rims are operated.)		Consult your dealer

Release the hoop	Blinking	Release the hoops during start-up. You can drive once the white driving profile indicator light comes on.
------------------	-----------------	---

Troubleshooting

Situations	Possible causes / solutions
System does not turn on with power button	<ul style="list-style-type: none"> Battery too discharged, connect iFlow battery charger. System in transport mode; see chapter on using the system box (iBox). System defect; consult your dealer, occupational therapist or advisor.
One of the control modules shows no signs of life while another is working.	<ul style="list-style-type: none"> Control panel not connected; check wiring.
Support does not drive with the desired behavior (not the desired acceleration, speed, braking or functionality such as hemiplegia mode).	<ul style="list-style-type: none"> The driving behavior and certain functionalities of the iFlow system are adjustable. Consult your dealer, occupational therapist, or advisor.
The system makes noise while driving.	<ul style="list-style-type: none"> Electric brakes are not properly adjusted; consult your dealer. Wiring touches the wheels; ensure the wiring is neatly concealed or consult your dealer, occupational therapist, or advisor.
A bump can be felt while driving.	<ul style="list-style-type: none"> By standing still in the same place for a long time, we have developed a flat side on wheels; moving with the wheelchair will improve the situation
System does not charge when connecting charger.	<ul style="list-style-type: none"> Make sure you are using an iFlow battery charger. Check that the power cord cable is properly connected to the iFlow battery charger. Clean the magnetic connector of the battery charger and control panel with a dry cloth. Check that the battery charger is connected to the wall socket
Horn not working.	<ul style="list-style-type: none"> The horn volume is adjustable; consult your dealer, occupational therapist, or advisor.

Maintenance matrix

The following matrix indicates which components should be checked regularly. If a problem occurs, contact your dealer to resolve the issue.

Element	Action	With every use	Weekly	Monthly	quarter	Annual
Battery charger	Check the functioning and condition of the cables	X				
Enclosures	Clean with a damp cloth		X			
Connections and cables	Check for damage and breaks		X			
Moving parts	Check function and clearance			X		
Battery	Fully charge			X		
Parking brake and freewheel levers or band brakes	Check functioning				X	
Complete system	Service					X

iFlow identification sticker

Wheelchair type identification sticker

Service and technical support

Consult your dealer, occupational therapist, or advisor for any questions, maintenance, or repairs. Always keep the following information at hand:

- The serial number of your iFlow system. You can find this on the System Box (iBox); see the chapter on using the System Box (iBox).
- Software version. You can find this in the status screen of the Joystick/Guide control.
- The wheelchair's serial number. You can find this number in the wheelchair's user manual.

Have the system checked and serviced by your dealer at least once a year. Preventive maintenance checks the proper functioning of all moving parts, as well as the electronics and the condition of the batteries. Repairs may only be carried out by your dealer or by Life & Mobility. Only genuine parts supplied by Life & Mobility may be used. Contact Life & Mobility for information about product safety notifications and recalls.

Warranty

Life & Mobility products are manufactured with care and are thoroughly inspected before leaving the factory. If a product does not meet your expectations, please contact the retailer from whom you purchased it. Life & Mobility complies with consumer protection laws and regulations in force in the country where the product was purchased.

For the warranty period, please refer to the Life&Mobility warranty overview. You can request this from our customer support department. The warranty does not apply to defects caused by failure to follow the user and maintenance instructions, abnormal use, wear and tear, damage, neglect, overloading, accidents caused by third parties, use of non-original parts, and defects not caused by the product.

This warranty replaces all other warranties, whether statutory or oral, except for written warranties from Life & Mobility.

Guarantees are only valid within the EU.

Afterword

This manual has been compiled with the utmost care, based on information known to Life & Mobility at the time of publication. The information in this manual pertains to the standard version of the product. Life & Mobility accepts no liability for damage that may occur as a result of changes or modifications to the product after it was first sold. Furthermore, Life & Mobility cannot be held liable for damage caused by work performed by third parties.

The information in this manual may not be reproduced for commercial purposes and/or published in any way without prior written permission from Life & Mobility.

In accordance with the Trade Name Act, the trade name and trademarks used by Life & Mobility are not considered to be free of charge.

No rights can be derived from this text. Subject to change.

© Life & Mobility – All rights reserved

Manufacturer Contact Information

Life & Mobility bv
Logistiekweg 7
7007 CJ Doetinchem
The Netherlands (NL)
T +31 (0)314 328 000

Postbus 304
7000 AH Doetinchem
The Netherlands
www.life-mobility.com
info@life-mobility.com

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