
MATCH

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Safety instructions	4
Introduction	8
Technical specifications	10
Product description	12
Settings and adjustments	16
Use	32
Maintenance	34
Transport & storage	35
Technical problems	39
Guarantee	41
Disposing of product & packaging	41

Safety instructions

General

Be careful:

- Never change the construction of the wheelchair.
- Use only original Life & Mobility parts and accessories. These are tested and safe.

Attention:

- Directions, safety symbols and instructions are attached to this wheelchair. They must not be covered over or removed. They must remain attached to the wheelchair throughout its lifespan and be clearly legible.

Allways notify Life & Mobility and the relevant authority in case of a serious incident.

Assembly and adjustment (see p. 16)

Warning:

- Before using the wheelchair, check that all its parts are properly mounted and secured.

Attention:

- Ensure that the push bars, the armrests, the lower leg length and the headrest, if applicable, are adjusted to the correct height.
- The foot supports must be at least 5 cm from the floor.
- In certain configurations, the wheelchair can be longer and/or wider than recommended according to the standard. This could make access to emergency exits more difficult. Use the adjustment options (push bar lower, seat angle upright, legrests and armrests at a different angle or remove them) in order to make the wheelchair more compact.

Safety instructions

Sitting in the wheelchair

Warning:

- Make sure that the occupant's fingers and limbs are not trapped between the frame and the wheels or between the wheels' spokes. Always be very careful when:
 - The wheelchair is tilted far backwards
 - The wheelchair is set to the highest seat level.
 - The occupant is sitting lopsided in the wheelchair.
 - The occupant has his/her arms hanging limp over the armrest.
 - A combination of the above occur.
- Prevent the feet of the user from ending up between the two foot supports. There is a heel strap or a central legrest available for this purpose.

Driving the wheelchair (See also p. 32).

Warning:

- The posture belt is not suitable to serve as a safety belt.

Attention:

- Pay attention to your hands when using the handrims if there is not much room at the side(s) of the wheelchair.
- Set the brake on the wheelchair before making a transfer.
- Use the brake only to block the wheels when standing still. The brake is not suitable for braking the wheelchair.
- The maximum pushing force required can differ for each configuration. Among other things, this depends on the weight of the occupant, the height of the push bars and the settings of the balance point.

Safety instructions

Transporting the wheelchair (See also p. 35).

Warning:

- In conformance with the covenant, the Match chassis was declared crash test safe in combination with specific moulded seats (ISO 7176-19:2008). In untested Match combinations, the chassis has not been declared crash test safe.
- Attach the chair in a taxi van using a 4 point securing system. Also use a 3-point safety belt to transport the wheelchair with the occupant.
- The posture belt is not suitable to serve as a safety belt.
- When transporting the wheelchair with the occupant, the wheelchair must be fitted with a headrest.
- Tilt the frame a maximum of 5degrees in a taxi van. This tilting angle does not apply to the moulded seat.

Crash test safe sticker

Slopes and other obstacles (See also p. 32).

Warning:

- Never drive on a slope at an angle exceeding 5 degrees.
- Do not drive on slopes if the wheelchair is tilted.
- When driving up a slope, the centre of gravity will shift and the wheelchair can tip backwards more easily.
- When driving down a slope, the centre of gravity will shift and the occupant can tip forward.
- Do not turn around on a steep slope. If the wheelchair is crosswise on a steep slope, the risk of tipping is significant.
- Always go down on steps backwards so the wheelchair cannot tip forward and to prevent

Safety instructions

the user from falling out of the chair.

- Assistants must never leave the wheelchair on a slope without applying the brakes.
- When moving the wheelchair by yourself, only use your hands on the handrims to brake and do not use the tyre brake. Otherwise, you risk losing control of the wheelchair.

Attention:

- When going down a slope, always keep your speed under control.
- The assistant must brake the wheelchair using the push bars.

Miscellaneous

Warning:

- Do not leave the wheelchair in direct sunlight: the upholstery may get too hot, which can burn the skin or make it hypersensitive.

Be careful:

- Do not allow the wheelchair to come in contact with seawater as seawater is corrosive and will rust the wheelchair.
- Do not allow the wheelchair to come in contact with sand, as it may penetrate into the moving parts of the wheelchair, causing unnecessarily fast wear and tear.
- Never use abrasive or aggressive cleaning products or solvents such as thinner, benzene, etc., to clean the wheelchair.
- Tightening the adjustment screw too far can damage the wheelchair.

Introduction

Dear user,

The Match is a tilting wheelchair developed particularly for use in combination with a moulded seat, for permanent wheelchair users who require a chair with (much) individual support and comfort.

The Match always provides the user and health care professional a customised solution, regardless of the dimensions or the type of seat support, the weight, body size and the health condition of the user.

The Match is resistant to spraying water and is suitable for indoor and outdoor use. Maximum occupant weight is 160 kg. This is stated on the CE sticker.

The Match received positive results tests conforming to the current requirements on European wheelchair standard EN 12182 and EN 12183, as well as with ISO 7176-19:2008 in combination with specific moulded seats as a crash test safe wheelchair (see page 31).

The Match therefore complies with the mandatorily applicable European REGULATION (EU) 2017/745 for Medical Devices.

In this user manual you will find everything you need to know about the Match, from safety instructions and control functions to maintenance tips. Sections important for your safety are marked in the text with a safety symbol.

CE sticker

Introduction

Warning:

- Not following the safety directions, instructions, warnings and maintenance recommendations may have a negative effect on the wheelchair performance and could result in physical injury or damage to the wheelchair or the environment.

Attention:

- Read this user manual carefully before using the wheelchair.
- Keep the user manual in a safe place for a possible future user.



We are happy to keep you informed on updates, user experiences and other useful information. You can read about the latest developments regularly on the www.life-mobility.com website.

Kind regards,

Management and staff of Life & Mobility

Technical specifications

Expected service life

Life & Mobility assigns a service life of 7 years to this product. This means that during this period, the product is fixable, safe and usable.

The statement above is only applicable if the user uses the product in the manner that Life & Mobility intended, marked with CE and including any original accessories.

Technical specifications

Abbrev.	Description	Measurements
OW	Moulded seat width Pushchair	42 - 67 cm (elbow 1 & 2)
	Moulded seat width Self-propelled	42 - 60 cm (at seat height column 1: max. moulded seat width 50 cm)
SH	Seat height*	30 - 60 cm (columns 1, 2, 3 & 4)
MD	Moulded seat depth	40 - 65 cm
WB	Wheel base	50 - 60 cm
WW rear	Wheelchair width Pushchair rear	55 - 70 cm (SW + 5 cm)
	Wheelchair width Pushchair front	60 - 75 cm (SW + 10 cm)
WW front	Wheelchair width front	60 - 75 cm (SW + 10 cm)
BP	Balance point	Adjustment length 15 cm
AH	Armrest height	17 - 40 cm
LW	Legrest width	
	Comfort legrest:	39 - 69 cm
	Comfort Solid legrest:	39 - 69 cm
BH	Back support tube height	52 cm
PH	Push bar height	80 - 150 cm (range adjustable, PH depends on SH)
LL	Lower leg length	
	Comfort legrest:	35-45 cm / 40 - 55 cm
	Comfort Solid legrest:	35-45 cm / 40 - 55 cm
	Central legrest:	depends on HMO
HMO	Height of moulded seat	Vertical distance top of assembly plate -bottom lower leg on the front of the moulded seat
-	Seat angle adjustment**	-5° -30° (range is adjustable)
-	Backrest angle adjustment	85° - 140° (range is adjustable)
-	Armrest depth	Adjustment length 8 cm
-	Total weight Pushchair	37,5 kg (not including moulded seat)
-	Total weight Self-propelled	39 kg (not including moulded seat)
-	Maximum occupant weight	160 kg
-	A CE-certified and crash test safe	EN 12182
	Match - moulded seat combination	EN 12183
	complies to	ISO 7176-19:2008
-	Ambient temperature	Life & Mobility recommends using the Match at an ambient temperature between -5°C and +40°C.

* **Attention:** At the lowest and highest seat heights, the wheel positioning relative to the shoulder joint is not optimal for driving of the wheelchair.

** **Attention:** Depending on height, moulded seat width, moulded seat shape and balance point, the adjustment range may be limited.

Product description

Checking the delivery

Before using the wheelchair, check that the delivery matches the product description (see p. 12-15). If you find any transport damage on delivery, please inform your dealer.

Options

In addition to the standard parts, the following options are available for the Match.

Legrests

- Central legrest
- Front bridge with Comfort legrests or Comfort Solid legrests

Armrests

- Armrest system (armrest elbow and T-piece in two sizes)
- Upholstered armrests
- Wide armrests
- P-shaped armrest (Self-propelled only)

Moulded seat


- Moulded seat-specific assembly system

Headrest

- Standard headrest
- Standard headrest small
- Snellius headrest
- Snellius headrest with Edison bracket

Brakes

- Drum brake system with brake pedals left and right
- Hub brakes with braking levers (standard on Self-propelled)
- Handbrake, operated by accompanying person

Product description

Worktop


- Plexiglas worktop
- Worktop made to size

Accessories


- Posture belt with assembly system
- Hoop covers
- Spoke protectors
- Luggage holder

Product description

Pushchair

Self-propelled

Product description

Product description

- A. Seat angle indicator
- B. Push bar
- C. Head support control lever
- D. Push bar height adjustment control lever
- E. Back support tube
- F. Back angle hinge
- G. Balance point settings
- H. Seat height adjustment column
- I. Rear wheel
- J. Taxi tie-down point left and right rear
- K. Anti-rollback support left rear
- L. Anti-tipping support left rear*
- M. Rear bridge
- N. Brake pedal right rear
- O. Front leg right front
- P. Taxi tie-down point left and right front
- Q. Front wheel
- R. Assembly plate
- S. Legrest (option)
- T. Armrest bend, T-piece and armrest (options)
- U. Moulded seat-specific assembly system
- V. Seat angle tilting control lever
- W. Seat frame
- X. Riding frame

Not pictured

- Seat angle limitation pin
- Key set
- Manual
- Other options & accessories

* With respect to preventive safety, the standard Match version is equipped with anti-tipping supports.

Settings and adjustments

The Match has numerous setting and adjustment options so the chair can be completely adjusted to the user's preferences. Life & Mobility recommends you consult your dealer for the settings. You can make the adjustments easily yourself, without tools.

SEAT FRAME

Moulded seat assembly

With this Match chassis, Life & Mobility - in collaboration with selected moulded seat suppliers - has made a crash test safe combination of a chassis and moulded seat possible, for the first time in history. Together with these selected vendors, we bear the responsibility for and guarantee crash test safe product combinations. This has been laid down in an agreement. Only a Match - moulded seat combination which complies with the agreement is crash test safe.

For Match - moulded seat combinations without an agreement, Life & Mobility passes the responsibility for the product combination to the offering party. The product combination will not have a crash test safe marking on it. In this case, take note of the following:

Attention:

1. A moulded seat must be mounted to the Match on the underside and back:
 - Mount the bottom plate of the moulded seat by placing four M8 bolts in the holes intended for this purpose in the Match mounting plate.
 - Attach the back section of the moulded seat to the Match's back tubing in four places using the mounting strips provided.
2. When attaching a moulded seat or upholstered parts, always use a bolted joint, never a threaded joint.

Settings and adjustments

3. Only use original Life & Mobility parts for leg, arm and head supports.
4. Never make machined modifications to Match wheelchairs, such as drilling, sawing or bending.
5. Always mount accessories using Life & Mobility mounting attachments (T-slot mount blocks) in the back tubing.

The Life & Mobility moulded seat assembly systems are specifically developed for various, selected moulded seats. Use the correct assembly system / the assembly system included in the delivery for the moulded seat you are going to use in combination with the Match chassis.

The assembly system consists of an moulded seat-specific hole pattern in the assembly plate (A) and an moulded seat-specific assembly system (B) on the back support tube (C).

Moulded seat-specific assembly system

The hole pattern intended for assembly of your moulded seat is marked with purple stickers.


Tip:

- Use the balance point settings for optimum accessibility of the assembly holes (see page 22).

Seat width

The seat width can be adjusted seamlessly for seat widths between 42 and 60 cm* by sliding the armrest elbow in or out.

Sticker marking
moulded seat hole pattern


* Max. moulded seat at seat height column 1 is 50 cm

Settings and adjustments

Adjust seat width

1. Loosen the socket screw (D).
2. Slide the armrest elbow in or out to the seat width desired.
3. Retighten the socket screw (D).
4. Repeat on other side.


Adjust seat width

Seat depth

The seat depth can be adjusted seamlessly for seat depths between 40 and 65 cm by sliding the legrest suspension.

Adjust seat depth

1. Loosen the socket screw (E).
2. Slide the legrest or the front bridge in or out until the desired seat depth is reached.
3. Retighten the socket screw.


Adjust seat depth

Seat height

A selection can be made from 4 different seat height-adjustment column lengths (F) to realise the desired seat height for the user.

Preferred seat height	Seat height column
30-38 cm	1
36-44 cm	2
44-52 cm	3
52-60 cm	4

Seat height column

Attention:

- At the lowest setting of seat height adjustment column 1 (the shortest one):
 - The clearance under the foot support can be limited by a small knee angle.
 - The seat angle tilt range forward is limited to 0° (instead of -5°).

Settings and adjustments


Tip:

- The seat height is still continuously adjustable for 8 cm in each seat height adjustment column. If the final seat height is still unknown, we therefore recommend selecting a longer seat height adjustment column.

Seat height adjustment


1. Loosen the 2 socket screws (A) in the frame.
2. Lift the seat frame from the riding frame with the seat height adjustment columns.
3. Loosen the socket screw (B) on the adjustment block in the seat height adjustment column (C).
4. Line up the bottom edge (D) of the adjustment block with the graduated marks for the correct seat height.

Seat height adjustment

Attention: At the adjustment block at the lowest seat height of the shortest column (column 1) a M5x10 set screw has to be used instead of a cylinder screw with washer.

5. Retighten the socket screw (8).
6. Repeat the above-mentioned for the other seat height adjustment column.
7. Lower the two seat height adjustment columns into the frame until the adjustment block rubs.
8. Secure the seating height adjustment columns by tightening the 2 socket bolts (A) of the two seat height adjustment columns again (first the lower socket bolts on both sides, then the upper ones).

Adjustment block

Seat tilt range

To utilise the full back tilt range of the Self-propelled Match, the following settings may be used. The tilt range is greater:

Settings and adjustments

- .. at a higher seat height (from SH 48 cm full tilt range).
- .. the further you adjust the balance point forward.
- .. if the track width matches moulded seat width.
- .. if the armrest depth is moved frontwards.
- .. when the P-shaped armrest is used.

Backrest angle adjustment:

The backrest angle of the seat frame is adjustable and can therefore be adapted to the backrest angle of the moulded seat.

Attention:

- After changing the backrest angle, the seat angle indicator has to be recalibrated.

Adjust the backrest angle:


1. Loosen the socket screws (E) on either side of the backrest hinge.

Attention: secure the back support tube!

2. Change the backrest angle.

Attention: the teeth of the two gear shift discs must completely mesh with each other.

3. Tighten the socket screws (E) on either side again.

Adjust the backrest angle:

Balance point settings

The seat frame can slide forward and backward from the riding frame, with the occupant and all seat support components in it. This allows for adjustment of the centre of gravity position of the user in the moulded seat with respect to the driving frame.

Settings and adjustments

Warning:

- Only adjust the balance point with the occupant in the moulded seat in horizontal position, thus 0° tilt.

Be careful:


- If the knee angle is very small, the legrest may rub against the damper suspension in the extreme rear balance point position.

Attention:

- When using the seat angle tilt, the position of the centre of gravity ends up further to the rear of the chassis, so you must therefore always make sure you follow our safety recommendations for balance point adjustment (the length of the arrow is not an indication for the amount of movement, only for the direction of movement).

Balance point setting recommendation on factory settings	-	+
Relatively low seat depth	←	
User physique: - leg weight relatively high - torso weight relatively high	←	→
Frontward and backward mobility relatively high		→
Back moulded seat relatively thick	←	
A relatively high number of slopes during use		→
Seat angle used most frequently relatively large		→

Balance point calibration


- Extra luggage influences the position of the balance point, keep this in mind when calibrating the balance point.


Settings and adjustments

Balance point settings

- Loosen the 2 socket screws (A).
Attention: Only unscrew the front socket screw entirely from the clamp block. The rear socket screw must never be **completely removed** as it limits backwards sliding of the assembly plate.
Tip: You can use the extension tool included in the delivery for the requisite Allen key to get good access to both socket screws of the balance point adjustment.
- Slide the seat frame to the front or the rear until the optimum position has been found.
Tip: Keep to the green or orange area of the graduation, which is visible through the peephole (B) on the side of the frame.
- Retighten the 2 socket screws (A).

Balance point calibration

Balance point calibration peephole


Seat angle

The seat angle can be adjusted easily by control lever (C) on the right push bar.

Attention:

- At the lowest setting of seat height adjustment column 1 (the shortest), the seat angle tilt range is limited to 0° forward.
- If the moulded seat is very wide, the seat angle tilt range can be limited due to the maximum track width of the rear wheels.

Adjusting the seat angle

Adjusting the seat angle

- Squeeze the right control lever (C) to disconnect the dampers.
- Set the seat angle required.
- Release the control lever again.

Settings and adjustments

Seat angle indicator

The user can indicate what the best seat angles or seat angle areas are with the seat angle indicator (D).


Attention:

- Make sure the seat angle indicator is properly calibrated before the optimum seat angle / seat angle areas are set.

Calibrating the seat angle indicator

- Tilt the wheelchair until the assembly plate is horizontal.
- Loosen the socket screws (A) on the bottom of the seat angle indicator a bit.
- Rotate the seat angle indicator (B) until the lower edge of the air bubble is on the 0° line.
- Retighten the socket screws.


Calibrating the seat angle indicator

Adjust the seat angle indicator

- Set the seat angle required.
- Use the tip of a pen to slide two of the four bolts (C) above and below the air bubble on the side of the housing
- Repeat this with the other 2 bolts for a second seat angle to be indicated.

Tip: If only 1 set of bolts is used, the other set can be slid completely upwards and downwards.


Adjust the seat angle indicator

Settings and adjustments

Seat angle range


The seat-tilt range is from -5 to 30°. This range can be limited on the front (limiting forward tilt) and on the rear (limiting backward tilt). To do so, move the 2 locking pins in the hole pattern on the discs under the seat frame. The seat angle can be secured in 5-degree increments by inserting the two pins in the hole directly above and below the lip/ridge (D).

Seat angle range hole pattern

Seat angle range hole pattern

- (E): Spare hole for 2nd pin. (This spare pin comes with the original user manual). This hole does not limit the seat angle.
- (F): Limitation max. 0° forward tilt (horizontal).
- (G): Limitation max. 5° forward tilt. The pin can always be left in this hole if complete range is desired.
- (H): Limitation max. 25° backward tilt.

Limitation max. 25° backwards

DRIVE FRAME

Recommended drive frame settings

When adjusting the wheel base (WB) take the moulded seat depth and agility of the user in forward/backward directions into account.

When adjusting the track width (TB) in the rear, take the moulded seat width, agility of the user in lateral direction as well as the manoeuvrability through doors, etc., into account.

When adjusting the track width in the front: take the ideal positioning of the feet of the user into account.

Wheel base

Settings and adjustments


Wheel base

The wheel base of the chassis is adjustable so a wheel base suitable to the user can be selected. This is defined by the seat depth of the seat support and the user's agility. A lower seat depth and more agility require a wider wheel base.


Wheel base adjustment

1. Loosen the socket screws (A).
2. Move the front leg with the wheel up to the graduated mark of the wheel base desired.
3. Retighten the socket screw with a torque of 14 Nm.
4. Repeat this with the other front leg; both front legs must be adjusted to the same wheel base using the graduated marks.

Attention: never pull the front legs beyond the red mark.

Track width

Wheel base adjustment

Track width

Rear wheels

The track width of the rear wheels can be adjusted by sliding the rear axle with the rear wheels in or out. The track width is defined by the seat width of the seat support and the user's agility.

Attention:


- When deciding the track width of the rear wheels, possible conflicts between moulded seat & the brake pedal or wheel must be taken into account for maximum seat angle tilting. In that case adjust the track width.

Settings and adjustments

Rear track width adjustment

1. Loosen the 4 socket screws (B).
2. Slide the rear wheel with the frame tube in or out to the track width desired.
3. Retighten the socket screws.
4. Repeat this with the other rear wheel; both rear wheels must be adjusted to the exact same width using the graduated marks.

Rear track width adjustment

Front wheels


The track width of the front wheels is adjustable by rotating the front legs inward or outward around the pivot point. This creates more room to position the feet between the front wheels.

Front track width adjustment

1. Loosen the socket screw (C).
2. Rotate the front leg with the front wheel to the track width desired using the setscrews (D) on the inside and outside of the hinge.

Attention: Turn the two set screws (D) on the inside and outside tightly against the rotating part of the leg; this is the rotation stop.


3. Retighten the socket screw (C).
4. Repeat this with the other front leg; both front legs must be adjusted to the same angle using the graduated marks (E).

Front track width adjustment

Wheel position

There are two bore positions for the two wheel sizes (see fig.). The seat height of the Self-propelled Match is here independent of wheel size.

Self-propelled

Settings and adjustments

SUPPORT


Armrests

The height, depth and width of the armrests can be adjusted (see 'seat width', pages 17 & 18).


Armrest height adjustment

You can choose from 2 heights of the armrest system in order to get the desired armrest height; it consists of an armrest elbow and a T-piece.

Preferred armrest height	Armrest system
17-25 cm	Armrest elbow 1 + T-piece 1
24-40 cm	Armrest elbow 2 + T-piece 2

Armrest height adjustment

Armrest depth adjustment

The height of the armrest can be seamlessly adjusted for each armrest system.

1. Loosen the socket screw (A) on the front of the armrest system a few turns.
2. Adjust the armrest to the correct height.
3. Retighten the socket screw.
4. Repeat this for the other armrest.
5. Both armrests can be adjusted to the same height using the graduated marks.

Armrest depth adjustment

1. Loosen the screws (B) entirely (other hole in armrest) or partially (using the slot).
2. Set the seat depth required.
3. Retighten the screws.
4. Repeat this for the other armrest.

Settings and adjustments

Legrests


The Match can be fitted with a central legrest or a front bridge to install Comfort legrests or Comfort Solid legrests.

Legrest depth adjustment

For central legrest depth adjustment see 'Seat depth', page 18.

The Comfort legrests and Comfort Solid legrests adapters can be adjusted independently and seamlessly.

1. Loosen the socket screw (C) on the bottom of the front bridge.
2. Set the preferred seat depth.
3. Tighten the socket screw.
4. Repeat this for the other side.


Legrest depth adjustment

Legrest height adjustment

Adjustment of the legrest height is decided by the height of moulded seat: the vertical distance top of installation plate - bottom lower leg on the front of the moulded seat.


Suspension (D) can be raised or lowered by rotating the front bridge or by installing a standard or raised adapter.

Legrest height adjustment

Legrest width adjustment

1. Loosen the socket screw (D) on the bottom of the front bridge.
2. Adjust to the desired width.
3. Tighten the socket screw.
4. Repeat this for the other side.

Legrest width adjustment


Settings and adjustments

Lower leg length adjustment

Lower leg length adjustment


1. Loosen the socket screw (F) a turn.
2. Set the preferred lower leg length.
Attention: take a distance of at least 5 cm to the floor into account.
3. Retighten the socket screw.


1

Securing the ankle hinge

1. Adjust the desired ankle angle.
2. Drill a hole in the hinge.
3. Secure the ankle hinge with a fastening pin.

2

3

Heel strap/continuous foot support

Attention:


- Prevent the feet of the user from ending up between the two foot supports. There is a heel strap or a continuous foot support available for this purpose.

Headrest

The headrest can be installed at exactly the right position by means of the standard hinge points.

Headrest height adjustment

1. Pull the control handle (A) up.
2. Adjust the headrest to the correct height.
3. Move the control handle (A) back down.

Headrest height adjustment

Adjust the headrest height to the individual position

After having been removed, the headrest can always be moved back to the right height immediately by adjusting the stop block.

1. Loosen the socket screw (B).
2. Slide the block (B).
3. Retighten the socket screw (B).

Settings and adjustments

Attention:

- The headrest control handle is fixed in the option where the headrest is not removable/adjustable.

The headrest can be easily removed to make transfers:

1. Pull the small, innermost handle (A) up.
2. Remove the headrest.

Push bar

Adjusting the push bar height

1. Pull the large, outermost handle (C) up.
2. Adjust the push bar to the correct height.
3. Push the handle against the back support tube again.

Adjusting the push bar height


Anti-tipping supports

The anti-tipping supports can be rotated by means of the pull knob (D). This can be used to realise an anti-tipping support position for outdoors (1), indoors (2) and completely switched off.

Attention:

- Make sure the pin is locked back in the hole again after adjustment.

Anti-tipping support adjustment

The anti-tip support - if properly adjusted - protects the wheelchair user from tipping backwards. The risk of tipping the wheelchair backwards is influenced by:

- the position of the user's center of gravity in the wheelchair (which is determined by factors such as axle position, backrest angle, seat angle, seat tilt, and weight distribution over the body (e.g., amputations)).

Settings and adjustments

- any restlessness or spasms of the user.
- the slope on which you are driving.

Adjust the anti-tip as necessary based on the above considerations.

Use

Use

Making transfers

The Match is a wheelchair for users with diverse handicaps. Transfers from the wheelchair to the bed, toilet or shower, for example, and vice versa, can be made in various ways. The lifting method depends on the degree to which the wheelchair user can cooperate. The Match enables both active and passive transfers, with the help of an assistant or lift.

Attention:

- Always set the brake on the wheelchair before making a transfer.

Transfer forward

1. Fold up the foot supports or turn the legrests.
2. Adjust the chair to a negative seat angle if necessary, as this facilitates standing up.


Tip:

- When transferring from a passive lift to the Match, first tilt the wheelchair with the occupant not in it so that the wheelchair user can be placed in the wheelchair in the best position.


Negotiating obstacles

1. Push the push bar down and push on the anti-rollback support with your foot at the same time.
2. Push the wheelchair onto the elevation (obstacle) and release the anti-rollback support again.
3. Push the wheelchair further and lift the rear wheels over the obstacle with the push bars.


Going up steps can be negotiated forward; riding down steps is better backwards.

Stairs

Steps

Gentle slope

Attention:

- Keep the speed under control when going down a gentle slope.

Brakes

The brakes of the Match are adjusted correctly as standard at the factory. The Match is equipped with drum brakes that can be operated by the foot pedals in the rear and on the side of the chassis. The foot pedals must be in down position when the chair is on the brake.

The braking force may have been reduced over the years or may no longer be the same for both wheels. In that case consult your dealer.

Self-propelled wheelchair


If you self-propel the Match, the brake is incorporated into the hub of the hoop wheel. Each wheel can be placed independently from the other on the brake either backwards or forwards by means of the braking lever.

Self-propelled brakes in neutral

The braking force may be reduced through intensive or prolonged use of the brakes. You then need to set the brakes.

Using an Allen wrench (3 mm), tighten the setting screw on the underside of the braking lever by a half turn.

Setting screw

Check that the brake power is now sufficient. Further adjust if necessary.

Maintenance

Minor maintenance (every 2 weeks)

You can perform minor maintenance on the Match yourself:

- Clean your wheelchair with a damp cloth, water and multi-purpose cleaner.
- Check the tyre pressure and tyre condition (wearing). Soft tyres make riding more difficult due to an increase in roll resistance.
- Use the valve nipple that comes with delivery to inflate the tyres (e.g. when using a bicycle pump). The proper tyre pressure is on the tyres.
- You must regularly spray the dampers and the damper blocking (the rectangular block at the end of the damper) with some grease (e.g. Teflon spray). This guarantees proper performance.
- Check that all fastenings are secure and retighten them if necessary. Bear in mind that, depending on the specific design of the orthosis (for example, in the case of split or modular orthoses), certain joints may require extra attention. Pay particular attention to critical transition points: areas where forces, movements or loads converge, such as the load-bearing connections between the seat plate and the back tube, and other transitions between different components or materials.
- Check the wheel tension on a regular basis (once per quarter, depending on the intensity of use) to make sure that your wheelchair continues to ride comfortably and safely. Ask your dealer for assistance if necessary.

Be careful:

- Never use abrasive or aggressive cleaning products or solvents such as thinner, benzene, etc., to clean the wheelchair.
- Tightening the adjustment screw too far can damage the wheelchair.

Major maintenance (annual)

Your dealer would be happy to handle the major maintenance for your wheelchair. That way you are certain you are getting the best possible use from your wheelchair. You can also go to your dealer for parts replacement.

Transport & storage


Illustration 1

Crash test safe sticker

Taxi van

The Match was designed to be secured in a forward direction in a taxi van. The Match was tested positively as a crash-test safe wheelchair in accordance with ISO 7176-19:2008. Thus the Match can be transported as a wheelchair with occupant in a taxi van (see also the indication sticker on the wheelchair).

However, the condition of the occupant and the doctor's recommendations are decisive for safe transportation. Life & Mobility does recommend the occupant make a transfer from the Match to a fixed seat if possible.

- The Match was dynamically tested in a forward direction, where the test dummy was secured in the taxi van with a hip belt and a shoulder belt. Its suitability during other collisions (such as a lateral collision) or other types of belts was not tested.
- The Match was tested in the various configurations. The critical configurations with Kelvin seating elements and Matrixx headrests are also 'crash-test safe'.
- The chair should be fitted with a Matrixx headrest.
- The Match should be attached to the floor of the taxi using a four-point fastening system in accordance with EN 10542-2.

Warning:

- In conformance with the covenant, the Match chassis was declared crash test safe in combination with specific moulded seats (ISO 7176-19:2008). In untested Match combinations, the undercarriage has not been declared impact test-safe. Therefore, it will not be provided with a crash test safe-sticker and taxi hook symbols.
- A wheelchair without a factory-installed crash test safe sticker is not a safely transportable wheelchair.

Transport & storage

- The positions of the four attachment points on the Match are indicated by a hook symbol sticker (see illustration 2).
- Easy access and manoeuvrability in motor vehicles can be significantly influenced by wheelchair size and agility. Smaller wheelchairs with a shorter turning circle generally offer better accessibility to the vehicle and manoeuvrability to a forward position. The Match was assessed positively for securability according to ISO 7176-19:2008, annex D.
- The Match should be fastened using anchoring systems that are suitable for the 4 attachment points.
- The wheelchair user must use a three-point safety belt in the taxi van. These safety belts must be attached to the taxi van.
- The hip belt must run low across the waist of the user. The angle between the belt and the floor surface must be between 30° and 75° (a larger angle is preferable) (see illustration 3).
- The hip belt must run obliquely across the front of the body where the hip and pelvis meet. The belt must fit the body well. The top belt must support the user entirely, as indicated in illustration 4.
- The diagonal belt must run over the shoulder and chest as shown in illustration 5.
- The belts must be tight enough but not uncomfortable for the user.
- The belt may not be twisted.
- The belt may not touch any wheelchair components (such as armrests or wheels) so the belt no longer supports the body of the user completely (see illustration 6).
- If possible, the Match must not be transported in tilted condition (max 5° backwards).
- The Match should be used with the user facing


Illustration 2


Illustration 3


Illustration 4


Illustration 5


Illustration 6

Transport & storage

forward. The wheelchair should be secured according to the securing system instructions from the manufacturer.

- The purpose of the three-point safety belt is to prevent head and chest injuries.
- To prevent other occupants from being injured, any worktop must be removed (if present). This worktop must be stored safely elsewhere in the vehicle.
- Accessories must be secured or removed if possible.
- After any accident, the manufacturer representative must evaluate whether the wheelchair can be used again.
- Changes to and/or replacement of the attachment points or structural parts and frame parts should not be implemented without first consulting the manufacturer.

Posture belt

Warning:

- The posture belt is not suitable to serve as a safety belt.

The posture belt is intended to position the occupant. Attach this by means of the Match posture belt in the back support tube.

Transport & storage

Warnings

- In conformance with the covenant, the Match chassis was declared crash test safe in combination with specific moulded seats (ISO 7176-19:2008). In untested Match combinations the undercarriage has not been declared impact test-safe.
- The Match has been dynamically tested in a forward direction, where the test dummy (= 103 kg) was secured with a hip belt and a shoulder belt.
- Secure the chair in a taxi van with a 4point fastening system. Also use a 3-point safety belt for transporting the wheelchair and occupant to reduce the risk of head and chest injuries.
- The sitting posture belt is unsuitable as a safety belt.
- Worktops should be removed and stored safely in the taxi.
- All accessories should be secured or removed and safely stored in the taxi.
- After every accident, the manufacturer's representative should examine the wheelchair to determine whether it can be used again.
- The fastening points or structural and frame components should never be changed and/or replaced without first consulting the manufacturer.
- Any batteries should be "spill-proof".
- The belt clasp must be positioned so that no wheelchair parts can come in contact with it during a collision.

Technical problems

Should your wheelchair be subject to direct inadvertent technical problems, we advise you to contact your dealer. The dealer will assess whether he can repair the chair or whether the wheelchair has to be returned to the manufacturer for repairs.

You may be able to solve a number of problems yourself.

The seat angle is no longer adjustable

- Check to see that the control cable is not kinked or stuck somewhere.
- Check to see that the damper moves when the handle is squeezed. If necessary, add some grease to the damper block.

The brakes do not have sufficient braking power or the brakes on either side do not have the same braking power

- You need to set the brakes, go to page 33.

Moving parts such as the push bar, height-adjustable armrest, legrest are secured

- The above-mentioned parts may have become stuck in the frame due to excessive forces on them. Carefully loosen the components using a plastic hammer.

Socket screws for the settings cannot be loosened

- The socket screws may have been tightened too far. Try to slip a pipe over the long end of the Allen key to increase the leverage. Now try to loosen the screw.
- The screw head may be stripped. Try to get hold of the screw with adjustable-joint pliers. If not successful, consult your dealer.

Technical problems

Attention:

- *All repairs other than those mentioned above should be done by a certified dealer.*
- *Sloppy repairs can endanger your health. Moreover, the warranty on the product will be cancelled.*
- *Use common sense and decide whether you yourself are capable of solving technical problems of the chair. Always consult your dealer if you have any doubt.*
- *When setting and replacing parts, you must comply with this manual and your dealer's advice.*
- *Rule out any hazardous situation to the assistant and the occupant of the chair!*

Warranty

Terms and Conditions of Warranty

The Life & Mobility products are made with care and are carefully inspected before leaving the factory. Should it emerge that a product does not meet expectations please contact the party where you purchased the product.

Life & Mobility respects the laws and regulations that protect the consumer and that apply in the country in which the product is purchased.

The warranty does not apply to defects caused by: failure to observe the instructions for use and maintenance, use other than normal use, wear and tear, negligence, overloading, third-party accidents, non-original parts used and defects not caused by the product.

This warranty replaces any other warranty, either stipulated by law or communicated verbally, except for anything guaranteed in writing by Life & Mobility.

Warranties are only valid within the EU.

Disposing of product & packaging

From an environmental perspective, we recommend that you dispose of the product and the packaging separately to as great an extent possible, allowing the used materials to be efficiently recycled.

In most cases you can turn in your used product to your wheelchair dealer. If this is not possible, you can turn in the push support at your local waste collection station.

Table des matières

This manual has been prepared with the utmost care and is based on information known to Life & Mobility at the time of publication.

The information in this manual applies to the standard version of the product. Life & Mobility does not accept responsibility for any damage that may occur due to changes or modifications to the product after selling it for the first time. Neither can Life & Mobility be held liable for any damage caused by third-party work.

The information in this manual may not be reproduced for commercial goals and/or made public by any means without prior written permission from Life & Mobility.

Pursuant to the Trade Names Act, the trade name and trademarks used by Life & Mobility may not be considered free.

No rights can be derived from this text. Subject to changes.

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© Life & Mobility – All rights reserved