NEOX

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Contents

This manual has been prepared with the greatest possible care and is based on information known to Life & Mobility at the time of publication.

Information for people with a visual impairment can be obtained from Life & Mobility. Contact details can be found on the last page.

The information in this manual applies to the standard version of the product. Life & Mobility accepts no liability for any damage that may occur as a result of changes or modifications to the product. Life & Mobility cannot be held liable for any damage caused by the work of third parties.

The information contained in this manual may not be reproduced for commercial purposes and/or made public in any way whatsoever without the prior written consent of Life & Mobility.

Under the Trade Name Act, the trade name trademarks used by Life & Mobility may not be considered free

No rights can be derived from this text. Subject to changes.

© Life & Mobility – All rights reserved

Foreword	4
Product Description	6
Newton seating system	7
Options	8
Technical Specifications	9
Safety Instructions	12
Settings & Customizations	17
Folding the wheelchair	30
Use	32
Transport	34
Maintenance	38
Technical Issues	41
Guarantee	42
Used wheelchairs/environment	43

Foreword

The NEOX is a multi-adjustable light weight foldable wheelchair with integrated Newton seating vision. The frame is adaptable to a wide range of body measurements, making it possible to precisely adjust it to the user. In combination with the right seating support of the Newton seating system, you can create a perfectly fitting wheelchair for every user.

The NEOX is suitable for both indoor and outdoor use. The maximum user weight is 136 kg. This is indicated on the CE sticker. The CE sticker is located on the left upper front frame tube under the seat mat.

The NEOX has been tested with a positive result in accordance with the applicable requirements of the European wheelchair standard NEN-EN 12183. And as a crash-tested wheelchair in accordance with the corresponding standards in ISO 7176-19:2022 (see page 34). The NEOX therefore meets the essential requirements of the European Regulation (EU) 2017/745 for Medical Devices.

Intended Use

The product is intended to be used to provide a mobility solution for people with walking difficulties.

Use

The product provides a seat support. In the seated position, the user can move around by propelling the product by manually turning the rear wheels, pushing off with the feet on the floor, or the product is pushed by an attendant. The product must be used by 1 person with a maximum user weight of 136 kg, in normal environmental situations.

The NEOX is a configurable product that can be configured according to the product overview of the product.

For the NEOX, the following scenarios can reasonably be described that do not meet the intended use, as described above:

- Failure to follow the warnings in the user manual
- Exceeding the weight and sizing limits as stated in the manual
- Failure to follow the instructions for safe use

Intended Users

The targeted users are people who have difficulty walking. Age ranging from 6 - 100+, weight up to 136 kg. The use of the wheelchair is not limited to people with a specific diagnosis.

Clues

Intended for users with limited to complete mobility limitations of the lower extremities.

Contraindications

Users with dementia or cognitive impairments should not sit down or get up from the wheelchair unassisted

Intended environment

The NEOX is intended for indoor and outdoor use on paved surfaces.

Important::

- Please read these instructions carefully before using the NEOX
- Failure to follow the safety regulations, instructions, warnings, and maintenance recommendations may adversely affect the operation of the NEOX and potentially result in bodily injury or damage to the wheelchair or the environment.

Note:

- Sections in this document that are relevant to your safety are marked with a safety symbol.

Product Description

Control of the delivery

Before you start using the NEOX, please check that the delivered wheelchair matches the product you ordered and the product description below. The product will be fully assembled and handed over to you by your local dealer, including all options ordered. Any personal adjustments can be made by your dealer. If you notice any damage during transit, please contact your dealer.

Life & Mobility trusts that your NEOX will meet your expectations and wishes you a lot of convenience from it.

Product description (standard)

1. Push handle
2. Newton FLX backrest
3. Newton FLX Contour Seat Cushion
4. T-piece with armrest and side shield
5. Tire brake
6. Footrest
7. Folding frame
8. Front wheel
9. Integrated taxi mount front and rear frame, left and right
10. Rear wheel with QR axle
11. Tilt support

Not pictured

- User manual
- Options/ Accessories

These instructions describe the standard version of NEOX. The NEOX can of course be expanded with many options and accessories, these are listed under 'Options'.

Newton seating system

A relaxed sitting posture is achieved through the combination of the appropriate Newton seat angles and individual body support. This gives the client more energy for daily activities and social interaction.

With a Life & Mobility wheelchair, the seat and back angle can be chosen to fit each user. In addition, the different seating support parts can be chosen modularly within the Newton seating system, depending on the user's impairment, activity and participation level. In this way, a suitable wheelchair can be customized for every user.

Advantages of Newton seating system

- 'Sitting in balance' due to the appropriate back angle and matching seat angle.
- Easy to adjust and "first-time-right" with the NEOX and Newton configurator and the step-by-step plan and instructional videos.
- High-quality back support through 3 adjustable support surfaces
- Optimal pressure distribution and seat stability thanks to 3D-shaped seat cushion.
- Interchangeability of parts between the Life & Mobility seats
- Modular product design with many options

Benefits for the user

- Optimal orientation through balanced head.
- Social interaction through horizontal line of sight.
- Minimal effort due to biomechanical balance.
- Freedom of movement of the arms due to improved body stability.
- No discomfort due to support of the individual back contour.
- No discomfort due to good pressure distribution and no frictional forces.
- Space for organs and breathing due to open sitting position.

More information and the latest news about the Newton seating system can be found on the Life & Mobility website www.life-mobility.com

Options

In addition to the standard components, the following (most common) options are available for the NEOX:

Different types of rear wheels

- Different types of tires
- Hoop covers
- Spoke Protectors
- One-handed drive
- Camber Axle Plate
- Auxiliary drive

Braking system

- Extended tire brake
- Push or pull brakes
- One-handed tire brake

Back angle

- Fixed back angle (foldable)
- Adjustable back angle

Newton seating system

- Newton FLX backrests
- Newton FLX Seat Cushions
- Newton Curve backrests
- Leg Support: Comfort Leg Rest, Footrest
- Armrest: PUR armrest, padded armrest, P-shaped armrest, swivel hemiplegia armrest
- Head Support: Small, Standard

Accessories

- Worksheets
- Anti-tipping supports
- Cane holder
- 2 or 4 point belts
-

For more information and the latest news about the different options, please consult the Life & Mobility Website www.life-mobility.com

Technical Specifications

Life & Mobility recommends using the NEOX at an ambient temperature between -5 °C and +40 °C.


Service and technical support


Consult your dealer if you have any questions about specific adjustments, maintenance, or repair. Always have the serial number handy. You can find this on the CE sticker on the wheelchair on the left upper frame tube under the seat mat.

Expected lifespan

Life & Mobility assigns a lifespan of 5 years to this product. This means that the product is repairable, safe, and usable during this period.

The above statement only applies if the user uses the product in the manner intended by Life & Mobility, with CE and including any original accessories

Specifications	
Seat width (SW)	34,5 – 57 cm
Seat depth (SD)	40 - 56 cm
Seat height (SH)	34 - 55 cm
Armrest height (AH)	20 - 30 cm
Lower Leg Length (LL)	38- 57 cm
Loaded seat angle	12° - 16°
Backrest height	43,5 - 62 cm
Backrest angle	-5° tot 30°
Weight (transport)	12,7 kg
Overall width (WW)	SW+18 cm
Folded average width (depending on selected seat width)	35,5 cm
Rear wheel	20" , 22" , 24" , 26"
Front wheel	5" , 6" , 7" , 8"
Camber of the rear wheel	0° / 3°
Brake type	Tire brake
Max. user weight	136 kg
Reflexes	ja
ISO7176-19 (crash test safe)	ja

Materials:

- Wheelchair Frame: Aluminium
- Seat and backrest: PE fabric, PU foam.
- Rear wheels: Aluminium, Fe PU
- Front wheel(s): Aluminium, Fe, PA6, PU
- Armrests: Aluminium, PA, PU
- Leg Rests: Aluminium, PA6
- Handlebar: Aluminium, EPDM

Corrosion protection: Powder coating and anodizing

Safety Instructions

Life & Mobility assumes no liability for damage or injury caused by negligence or non-compliance with the safety instructions.

Depending on the specific conditions or accessories used, additional safety instructions may be required.

Always inform Life & Mobility and the competent authority in the event of a serious incident.

General

- Never change the construction of the wheelchair.
- Use only original Life & Mobility parts and accessories. These have been tested and are safe.
- Do not add extra weight to the wheelchair.

Markings on the wheelchair

- Markings, safety symbols, and instructions are displayed on this wheelchair. They should not be covered or removed. They must remain on the chair and be clearly visible throughout the life of the wheelchair.

Assembly and adjustment

- Before using the wheelchair, check that all parts of the wheelchair are properly assembled and secured.
- Make sure that the push handles, the back height, the armrests, the lower leg length and possibly the headrest are adjusted to the correct height.
- It is advised to maintain a minimum distance of 5 cm between the leg rests and the floor.
- Check the dimensions and adjust them optimally as described in this manual.

Sitting in the wheelchair

Do not allow the occupant's fingers and limbs to be trapped. Be particularly vigilant when:

- The occupant is skewed in the wheelchair.
- The occupant does not have full control over his/her arms and/or legs.

Use of the wheelchair

(see also page 32 for more information)

Warning:

- *The NEOX is not equipped with an anti-tipping support as standard. Check that tilting the wheelchair backwards does not suddenly allow it to tip backwards. In that case always add the anti-tippers for additional safety.*
- *When driving up a slope, the center of gravity will move backwards and the wheelchair can tilt backwards more easily.*
- *When going down a slope, the center of gravity will move forward and the user can fall forward from the wheelchair more easily.*
- *As a wheelchair user, you are vulnerable in traffic. Make sure that you are always seen and noticed by other road users. Obey the applicable traffic rules.*
- *Avoid lonely routes, make sure that help can be provided in case of an emergency.*
- The occupant or attendant must be alert to the possibility of the wheelchair tipping backwards in the following cases:
 - - When the balance position is set to active.
 - When the angle of the backrest is large
 - When driving up a slope.
 - When someone is standing on the anti-tip aid support and/or on the anti-tipping support.
- When self-driving the wheelchair, pay attention to your hands if there is not much space on the side(s) or brakes of the wheelchair.


- Keep in mind that there are fixed and loose parts sticking out, in and on the wheelchair. These can cause damage to the environment and to the wheelchair.
- Companions must never lose sight of the wheelchair on a slope without applying the brake.
- Never take obstacles higher than 10 cm alone.
- Apply the wheelchair brake before getting out of the chair or transferring.
- Only use the brake to lock the stationary wheels. The brake is not intended to slow down the wheelchair.
- When braking, heat is generated by the friction between your hands and the hoops. If you brake too hard, your hands can get too hot.

Transportation of the wheelchair without an occupant (see also page 34 for more information)
The NEOX is equipped with a folding frame and is therefore easy to transport in a car. Keep the following points in mind:

- Be careful not to pinch your fingers when folding the wheelchair.
- Always carry the wheelchair by the frame, never by the leg rest, armrests, wheels or the designated handle on the upholstery.
- Make sure the wheelchair is fully unfolded before using it.
- Do not place the wheelchair in direct sunlight. The upholstery may become too hot and cause burns or skin irritation.

Transporting the wheelchair with the occupant in it in a taxi bus (see also page 34 onwards for more information)

The NEOX is designed to be secured in a taxi van in a forward direction. The NEOX has tested positive as a crash test-safe wheelchair in accordance with ISO 7176-19:2022. So it can be transported with the

occupant in it in a taxi bus. See also the sticker on the wheelchair that is located on the upper frame tube on the front right under the seat mat.

- Attach the wheelchair with occupant to a taxi bus with an approved 4-point restraint system.
- In addition, use an approved 3-point safety belt for the safe transport of the wheelchair occupant.

Note:

- The normal sitting position belt is not intended to be used as a seat belt.

Ramps and other obstacles

(see also page 32 onwards for more information)

- The wheelchair is more likely to tip over on a slope than on a flat surface.
- Always move backwards off a sidewalk to prevent the wheelchair from tipping forward and the occupant from falling out. Never go down the stairs alone in the wheelchair.
- Never drive the wheelchair up a slope with an angle of more than 7°.
- Ask someone to help you if you're not sure if you can reach the top of the slope by your own strength.
- Lean or push your body toward the top of the slope. This shifts the center of gravity of the wheelchair forwards or backwards and reduces the risk of the chair tilting forward or backward.
- The person pushing the wheelchair must continue to push the chair to prevent it from starting to move backwards under its own power.
- Do not turn on a steep slope. When the wheelchair is perpendicular to a steep slope, the risk of tipping over is high.
- Keep your speed under control when descending slopes.


Settings & Customizations.

- The attendant must reduce the speed of the wheelchair using the push handles.
- If you are using the wheelchair yourself, slow down the wheelchair with your hands only on the hoops. Never use the brake, as you may lose control of the wheelchair.

Others

- Drive carefully on roads that are slippery due to rain, ice or snow!
- Do not allow the wheelchair to come into contact with seawater: seawater is aggressive and will affect the wheelchair.
- Do not allow the wheelchair to come into contact with sand: sand can penetrate into the rotating parts of the wheelchair, causing unnecessarily rapid wear and tear.
- Do not place the wheelchair in direct sunlight: the upholstery may become too hot and cause sunburn or skin sensitization.
- If you are under the influence of substances that may affect your ability to drive, you should not drive the wheelchair.
- You must have sufficient eyesight to be able to drive the wheelchair safely in the respective use situation.
- Be careful not to let any items of clothing hang loose. These could get between the wheels.
- Never lift the wheelchair by the armrests, as this can lead to dangerous situations.


The NEOX has a wide range of settings and adjustment options that can be used to adapt the wheelchair to all user requirements. The NEOX instructional videos that you can find online on our website or directly via the QR code on the right, explain all the settings clearly. The NEOX configurator provides support to achieve the right individual seat settings for your client. Life & Mobility recommends that you consult your dealer regarding the correct settings. You can make the simple adjustments yourself, without the use of tools.

Frame height

Use the QR code above or check out our website for the complete NEOX frame height chart with the right setup advice for your client.

Based on the lower leg length, way of locomotion, activity level, user profile and desired front and rear wheel size, an exact recommendation is given for the positioning of the front and rear wheel.

Frame height is measured at the front of the wheelchair from the top of the frame's seat tube to the ground.

Frame height front wheelchair with telescope

1. Loosen the 2 Allen screws (A).
2. Slide the front wheel telescope in or out of the front module until the front of the wheelchair has the desired frame height.
3. Place the nut plate and Allen screws (A) in the correct position and tighten
4. Repeat for the other front wheel.

Frame height front wheelchair without telescope:

1. Remove the 2 Allen screws (A) and the front wheel.
2. Place the front wheel at the correct height on the front module to achieve the desired frame height.

- 3. Retighten the 2 Allen screws (A).
- 4. Repeat for the other front wheel.

Adjusting the height position of the front wheel in the fork


- 1. Loosen the front wheel axle nut and remove the axle.
- 2. Place the front wheel back in the desired height position according to the NEOX frame height table or the advice of the NEOX configurator.
- 3. Replace the shaft and tighten the nut.

Positions 1 to 5 correspond to the height of the wheel in the clevis.

Frame height at the back of the wheelchair

- 1. Remove the rear wheel by pressing the button in the middle of the wheel.
- 2. Loosen the 2 axle plate bolts (A & B) and remove the axle plate from the rear module.
- 3. Place the axle plate in the correct height position on the rear module according to the NEOX frame height chart or the advice from the NEOX configurator, insert the 2 Allen screws (A & B) and tighten them. (tightening torque 10Nm)
- 4. Insert the rear wheel back into the axle bushing (C).
- 5. Repeat for the other rear wheel.

Note:

- Make sure the left and right axle plates are mounted at the same height.
- After adjusting the frame height, always check the head angle adjustment and the adjustment of the tire brakes (see page 23 onwards)

Seat width and back width


A. Frame Settings

The seat width of the wheelchair is adjustable by the width-adjustable folding frame and the width-adjustable upholstery. And can be adjusted in steps of 2.5 cm from 34.5 cm to 57 cm.

The size 1 folding frame has an adjustment range of 34.5 to 49.5 cm. And the folding frame size 2 has an adjustment range of 42 to 57 cm.

Adjusting the seat width of the folding frame

- 1. Remove the seat cushion.
- 2. Remove the cover caps (X) and the upholstery screws on the front and rear of the seat frame tubes and the top and bottom of the back tubes
- 3. Slide the upholstery rails with the upholstery forwards out of the seat frame tubes and up out of the back frame tubes.
- 4. Loosen the 6 Allen screws (A to F) of the folding frame.
- 5. Extend the folding frame to the desired seat width (teeth for every 1.25 cm, resulting in increments of 2.5 cm at a time). See the table on the right for which number of the folding frame should be visible for which seat width.

1	W34,5
2	W37
3	W39,5
4	W42
5	W44,5
6	W47
7	W49,5
8	W52
9	W54,5
10	W57

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


- Make sure that the left and right sides of the folding frame are set to the same number of teeth.
- 6. Retighten the 6 Allen screws (At/m F). (tightening torque 10 Nm).

B. Seat upholstery adjustment


1. Slide the two upholstery rails (A) into the correct seat width size indication (C) on the end of the seat width tunnels (B) of the upholstery.
2. Fold the end of the upholstery around the rails. Make sure that the rails themselves are not folded over, so the 'fold line' is outside the rails, see also the size indication numbers on the upholstery (C)
3. Place the left and right upholstery tools on the end of the frame tubes.
4. Fold the folding frame slightly, just like the back frame, and slide both rails (with the upholstery) into the tubes of the seat frame at the same time, using the upholstery tool
5. Screw the upholstery screws on the front of the upholstery into the tunnels where the upholstery rails are located.
6. Replace the cover caps (X) at the ends of the seat frame tubes.

**C. Back upholstery adjustment**

1. Slide the upholstery tendon (A) into the correct back width size indication (C) on the back width tunnels (B) of the upholstery.
2. Fold the end of the upholstery around the rails. Make sure that the rails themselves are not folded over, so the 'fold line' is outside the rails, see also the size indication numbers on the upholstery (C)
3. Fold the folding frame slightly, just like the back frame, and slide both rails into the tubes of the back frame at the same time, using the upholstery tool
4. Tighten the upholstery screws on both sides of the front of the upholstery in the tunnels where the upholstery rails are located.
5. Replace the cover caps (Y) on the top and bottom of the left and right back tube.

**Seat depth**

The seat depth can be adjusted by stepless sliding the rear module over the side frame.

1. Loosen the 4 Allen screws (A to D) on both sides.
2. Slide the rear module over the side frame to the desired seat depth position (see seat depth indication on the inside of the upper seat frame tube).
3. Retighten the 4 Allen screws (A to D) on both sides. (tightening torque 10Nm)

**Warning:**


- Make sure the rear module is mounted at the same depth on the left and right.
- If the position of the back frame prevents you from sliding the rear module to the correct position, adjust the back depth position according to the instructions on p. 22.

Seat depth adjustment range per side frame:

- Short frame: SD 40 cm - SD 48 cm
- Long frame: SD 46 cm - SD 56 cm

Back Depth Position / Fine Adjustment of Seat Depth

The position of the back support can be adjusted extra precisely by moving the back tubes over the seat frame tube (inside the rear module). The middle position between the rear module is the standard position.

1. Loosen the 2 Allen screws (C) of the back depth adjustment.
2. Slide the back tubes to the desired depth position.
3. Retighten the 2 Allen screws (C). (tightening torque 10 Nm)

**Warning:**

- Make sure the back tubes are mounted at the same depth on the left and right.

Adjustment of the back angle (in case of an adjustable back angle)

- Loosen the Allen screw (A) on both back corner hinges on the back.
- The indicator line on the side of the back angle hinge shows the standard Newton hip angle of 103 degrees (angle between backrest (LT plane) and seat).
- Retighten both Allen screws (A).

Warning:

- The backrest angle adjustment may only be carried out without an occupant.
- Make sure the left and right back tubes are mounted at the same angle.

Backrest Height Adjustment

- Loosen the 2 Allen screws (B) from the height position of both back tubes.
- Slide the back tubes to the correct LT plane height (see indication from 0 to 8 on back tubes) The NEOX configurator indicates the correct height position of the backrest based on your client's LT flat height.
- Retighten the 2 Allen screws (B) on both back tubes. (tightening torque 10 Nm).


Note:

- Make sure the left and right back tubes are positioned at the same height.
- Check that the back tube is properly fastened before anyone uses the chair.

Balance Point

The balance point of the NEOX can be adjusted to your client's needs. The NEOX configurator gives you advice on a suitable balance point setting for your client.

- Loosen the large axle bushing nut (A) on the inside of the axle plate and place the axle bushing in the correct hole position.
- Retighten the nut on the inside of the axle plate.


Warning:


- Make sure the axle bushing is in the same axle hole position on both sides of the wheelchair.
- Never use the front axle hole position for seat angle B.
- The front axle hole position results in the most active, light-moving wheelchair, but also the most tilt-sensitive wheelchair.
- The wheelchair is least sensitive to tipping in the rear axle hole position.
- When changing the seat height, seat angle and axle plate position, the tire brake and head angle of the front wheel must be readjusted if necessary, so that the handling of the seat is optimal again and the braking effect is sufficient.

Adjusting the Head Angle

If the rear wheel position is moved up or down with the axle plate, or the front wheel is moved up or down, the angle of the front fork changes relative to the floor. In order to maintain the favourable, agile handling characteristics of the NEOX, the castor axle must be perpendicular to the floor.

The axis around which the castors rotate, the castor axis, must always be placed perpendicular to the floor (this is at an angle of 90 degrees).

- Loosen the 2 Allen screws (A) slightly, so that the ball head (B) can just be rotated by hand.
- Then place a spirit level or a smart phone with a level app against the vertical side of the head-set. (B)
- Now insert an Allen wrench into the hole of the eccentric (C) and twist around the Allen screw. The angle of the fork will change.

- Adjust the angle of the fork relative to the floor to 90 degrees.
- Retighten the Allen screws..

**Warning:**


- As this is difficult to see, you can check this by driving the chair forward and backward slightly, with the caster fully swinging from back to front. If the castor axle is vertical, there is no measurable or visible difference in height between the frame tubes and the substrate during this operation.

Adjusting tire brakes

When changing the seat height, seat depth, seat angle and balance point, it is necessary to readjust the tire brakes.

Rough adjustment


- Loosen the Allen bolt (A) on the underside of the brake clamp.
- Slide the brake forward or backward through the seat frame tube until the desired position in relation to the tire is reached.
- Retighten the Allen bolt (A) on the underside of the brake clamp.

**Warning:**

In certain configurations, the tire brakes are mounted on the front module. Then it is not possible to use the rough adjustment, only the fine adjustment.

Fine adjustment

- Loosen the Allen bolt (B) at the top of the brake clamp.
- Set the band brake to the desired position.
- Retighten the Allen bolt (B) at the top of the brake clamp.

**Adjusting the height of a push handle**

- Unlock the quick-release lever (A) by moving the lever upwards
- Adjust the push handle to the desired height.
- Lock the push handle by placing the quick-release lever (A) vertically against the frame

Armrests

The armrests are adjustable in height, depth and width.

Inserting the armrests

Insert the plug-in tube with the armrest straight down into the side holder.

**Note:**

- Make sure the side shield is on the seat side, not the wheel side.

*Removing the armrest*

Lift the armrest (with side shield) vertically up from the side holder.

Adjusting the height of the armrest

- Loosen the Allen screw (A) on the front of the armrest.
- Adjust the armrest to the desired height.
- Retighten the Allen screw (A).

Adjusting the depth position of the armrest

- Loosen the Allen screw (B).
- Set the desired depth.
- Retighten the Allen screw (B).

Adjusting the width position of the side shield

1. Lift the T-piece armrest out of the side holder.
2. Remove the height stop from the T-piece by loosening the Allen screw.
3. Slide the mounting blocks (A) out of the insertion tube of the armrest and split it into 4 loose blocks.
4. Slide the 2 threaded mounting blocks into the plug-in tube.
5. Place the 2 additional blocks on the 2 threaded mounting blocks and place the side shield on top.
6. Assemble this assembly with the 2 bolts supplied.

Fix the armrest

1. Slide a thin square nut down behind the height stop.
2. Insert a countersunk bolt (D) into the hole and tighten the bolt into the square nut.

Adjusting the depth of the side shields.


1. Remove the 2 countersunk bolts (C).
2. Place the side shield in the desired position.
3. Retighten the 2 countersunk bolts (C).

Adjusting the height of the side shields.

1. Loosen the 2 countersunk bolts (C).
2. Slide the side shield up or down.
3. Retighten the 2 countersunk bolts (C).

Replace the armrest

1. Remove the 2 countersunk bolts (E) on the bottom of the armrest.
2. Remove the armrest and replace it with a new one.
3. Retighten the 2 countersunk bolts (E).

Leg rests

Swivel and/or remove the leg rests

1. Pull the lever (C) upwards and rotate the legrest outwards or inwards at the same time.
2. If the leg rest is turned outwards, it can be lifted straight up from the frame.

Positioning the leg rests

1. Rotate the leg rest a quarter turn outwards from the mounted position and lower it into the opening provided.
2. Rotate the leg rest forward until it clicks into place.

Warning:

- If a push brake is fitted, make sure that the leg rest cannot unintentionally release the tire brake when swinging outwards.

Adjusting the lower leg Length

1. Loosen the Allen screw (D)
2. Slide the insert tube further in or out of the suspension tube until the desired lower leg length is reached.

Note:

- It is advised to keep at least 5 cm of vertical distance between the footplate and the ground.

3. Tighten the Allen screw (D).

Adjusting the Knee Angle

1. Remove the countersunk bolt (A) in the ankle joint and footplate
2. Loosen the Allen screw on the underside of the leg rest (B).
3. Remove the single pivot point from the plug-in tube and rotate it 180 degrees.
4. Replace the single pivot point in the plug-in tube and re-tighten the Allen screw (B).
5. Replace the footplate against the pivot point and retighten the countersunk bolt (A) in the ankle joint.

**Note:**

- Check that there is no conflict between the front wheels and the leg rests.

Anti-tipping supports

The anti-tipping supports (optional) can be temporarily swivelled forward to make it easier, for example, to climb a curb or a high threshold.

1. Push down the horizontally spring loaded part of the anti-tipping support (A).
2. Turn the anti-tipping support inwards and forwards.
3. Repeat for the other side.

**Warning:**

- Inform the occupant when the anti-tipping supports are folded away and ensure that the anti-tipping supports are pushed back to their original position after the obstacle has been cleared.


Adjusting the height of the anti-tipping support:
Rough adjustment

1. Loosen the 2 Allen screws that secure the anti-tipping support to the rear module.
2. Move the anti-tipping support to the other mounting position
3. Retighten the 2 Allen screws (C).

Fine adjustment

1. Remove the Allen screw (B) and the bushing.
2. Adjust the anti-tipping supports to the desired height.
3. Replace the dowel and tighten the Allen screw (B).

The anti-tip support - if properly adjusted - protects the wheelchair user from tipping backwards. The risk of tipping the wheelchair backwards is influenced by:


- the position of the user's center of gravity in the wheelchair (which is determined by factors such as axle position, backrest angle, seat angle, seat tilt, and weight distribution over the body (e.g., amputations).
- any restlessness or spasms of the user.
- the slope on which you are driving.

Adjust the anti-tip as necessary based on the above considerations.

Adjustment of the camber of the rear wheel

Optionally, an axle plate with 3 degrees of camber is available. To make use of this option, follow the steps below.

1. remove the wheel (QR) and axle bushing (C) from the frame
2. Remove the wheel cover by loosening the 2 Allen screws (A and B).
3. Insert the camber wheel plate and secure it with the two Allen screws.
4. Put the axle bushing back in the desired position.
5. Insert the wheel

Check that there is enough space between the wheel and the side of the frame.

Folding the wheelchair

Folding the wheelchair

Remove the seat cushion

Remove the seat cushion from the seat mat.

Remove the leg rests (optional)

Remove the leg rests according to the instructions on pg. 27.

Remove the armrests with side shields (optional)

Remove the armrests and according to the instructions on p. 25.

Remove rear wheels (optional)

Press the button in the middle of the wheel (the QR axis) while pulling the wheel out of the frame.

Folding push handle (optional)

Unlock the push handles with the quick release (A) and turn the push handles downwards.

Unlocking the tension rod

Unlock the folding tension rod (B) by pulling up on the middle part, which will cause the rod to fold. Or unlock the non-folding tension rod on one side of the rod so that it hangs down.

Fold

The chair is now ready to be folded

1. Push the green foot pedal (C) up. The folding frame will begin to fold.
2. Pull up the upholstery with the fabric handle.
3. The wheelchair is now getting smaller.

Fold out

1. Push the two sides of the wheelchair apart.
2. Push the green foot pedal down.
3. Lock the tension bar and put all parts back on the wheelchair

Warning

- Be careful not to get your fingers caught between the frame.
- Make sure the seat is fully unfolded before riding it.
- Check that all parts are secure before you or your client return to your wheelchair.

Use

The NEOX is a wheelchair for users with different disabilities and needs. Transfers from the wheelchair to, for example, a bed, toilet or shower and vice versa can be made in different ways. The extent to which the wheelchair user can cooperate determines how he or she should be lifted. Both active and passive transfers are possible with the NEOX, possibly with the help of the caregiver or with the help of a lift.

Getting in and out of the wheelchair yourself

Warning:

- Apply the wheelchair brake before entering or exiting.
- Do not use the brake as an aid to stand up.
- Flip up the footplates.
- Operate the leg rest pull lever and rotate the leg rest outward and/or remove the leg rest.

In the case of a lateral transfer, you can remove the sides if necessary.


Overcoming hurdles

When overcoming obstacles independently, Life & Mobility advises extensive practice under the guidance of a therapist.


Warning:

- When going up a slope, the center of gravity will shift backwards, making it easier for the wheelchair to tilt backwards.

Life & Mobility advises less experienced users to adjust the wheelchair less actively and to equip it with anti-tipping supports.

Figuur 1

Figuur 2

Overcoming hurdles with a companion:

- Push down with your foot on the tilt support and push the push handles down, so that the front of the wheelchair lifts off the ground.
- Push the wheelchair onto the elevation (obstacle) and release the tilt support.
- Push the wheelchair further and lift the rear wheels over the obstacle using the push handles.

Driving up a sidewalk can be done forward, driving down a sidewalk is better done backwards. (Figure 1).

Please note:

- Keep your speed under control on slopes. (Figure 2).

Brakes

The brakes of the NEOX are set up correctly as standard.

- The correct way of braking is by pulling both the left and right brake levers backwards (in the case of pull control) or pushing forward (in the case of push control).
- The brake is disengaged by applying a slight forward (or backward) pressure to the brake lever. The brake is brought back to the correct position by a spring.
- The brakes are intended as a parking brake only, not to slow down the wheelchair.

Transport

Taxi bus

The NEOX is designed to be secured in a taxi bus in a forward direction. The NEOX has tested positive as a crash test-safe wheelchair in accordance with ISO 7176-19:2022. The NEOX can therefore be transported as a seat with occupant in a taxi bus (see the sticker on the wheelchair on the upper frame tube, front right under the seat mat).

However, whether the occupant of the wheelchair can be transported safely depends on his or her state of health. The user's physician should be consulted. Life & Mobility recommends moving the occupant from the NEOX to a fixed seat if possible. The wheelchair must then be securely secured in the vehicle.


- The NEOX was tested, where the test dummy was secured with a 3-point belt. The suitability for collisions or other types of belts than those mentioned in the above standard has not been tested.
- The NEOX has been tested in different configurations.
- The NEOX must be equipped with an approved* 4-point hold-open system in accordance with EN 10542-2, attached to the floor of the minibus.

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- If there is no crash test safe sticker on the wheelchair from the factory, this is not a wheelchair that can be transported safely.
- Only lift the wheelchair by fixed frame parts and not by removable accessories such as armrests, leg rests, seat, backrest, etc.
- The positions of the four attachment points on the NEOX are indicated by a hook symbol sticker (see Figure 1).

These points must be used to secure the wheelchair.

hook symbol back

hook symbol front


Figure 1


Figure 2


Figure 3


Figure 4

Easy access and maneuverability in motor vehicles can be significantly affected by the size of the wheelchair. Smaller wheelchairs with a shorter turning radius generally offer better vehicle accessibility and better maneuverability to a forward position. The NEOX has a good rating with regards to fixability according to ISO 7176-19:2008 annex D.

- The Neox must be transported with the user in a forward position. The wheelchair must be secured in accordance with the instructions of the restraint manufacturer.
- The NEOX must be secured with anchoring systems that are suitable for the 4 attachment points.
- The wheelchair user must use an approved 3-point seat belt in the taxi bus. These seat belts must be attached to the taxi bus.
- The lap belt should run low over the user's waist. The angle between the belt and the floor area should be between 30° and 75° (a larger angle is more desirable) (see Figure 2).
- The lap belt should run across the front of the body at the point where the hip and pelvis meet. The belt must fit snugly to the body.
- The upper belt must fully support the user as shown in Figure 3.
- The diagonal belt should run over the shoulder and chest as shown in Figure 4.
- The belts must be sufficiently tight and not uncomfortable for the user
- The belt must not be twisted.
- The belt must not lie against wheelchair components (such as armrests or wheels), which means that the belt no longer fully supports the user's body (see Figure 5).
- The approved 3-point safety belt is intended to prevent injury to the head and chest.

- To prevent injury to other occupants, a worktop (if any) must be removed. This worktop must be safely stored elsewhere in the vehicle.
- If possible, the accessories should be secured or removed.
- After any accident, the manufacturer's representative must assess the wheelchair to see if it can be used again.
- Modifications and/or replacements to the tie-down points or structural parts and frame parts should not be carried out without prior consultation with the manufacturer.

Sitting position belt

The sitting position belt is intended to permanently position the occupant in the desired sitting position, e.g. in the event of restlessness, by difficulties maintaining a correct position, or when overcoming an obstacle. Attach it to either side the wheelchair frame using the associated mounting brackets (option).

(with exclamation mark symbol)

Please note:

- The seat position belt is unsuitable as a seat belt.

Car

In addition to a taxi bus, the NEOX can also be transported folded in the trunk of a passenger car.

See page 27 for instructions on how to fold the wheelchair.


Figure 5

Warnings

- In accordance with the covenant, the NEOX chassis in combination with specific orthoses has been declared crash test safe (ISO 7176-19:2008) and provided with CE marking.
- The NEOX was dynamically tested in the forward direction, where the test dummy (= 76 kg) was secured with a lap belt and a shoulder belt.
- Attach the seat in a taxi bus with a 4-point restraint system. Also use a 3-point seat belt to transport the wheelchair with occupant to reduce the risk of injuries to the head and chest
- The seat position belt is unsuitable as a seat belt.
- Worktops must be removed and stored safely in the taxi.
- All accessories must be secured or removed and safely stored in the taxi.
- After any accident, the manufacturer's representative must assess the wheelchair to see if it can be used again.
- Modifications and/or replacements to the tie-down points or structural parts and frame parts should not be carried out without prior consultation with the manufacturer.
- Any batteries must be "spill-proof".
- The buckle of the belt must be positioned in such a way that wheelchair parts will not make contact with it during a collision.

Maintenance

Maintenance

You can do minor maintenance on the NEOX yourself.

Cleaning and disinfection

Cleaning the back cushions and seat cushions can be done by washing them according to the washing instructions on the product. The frame, push handles and armrests can be cleaned with a damp cloth, water and commercially available, non-abrasive detergent with a pH between 5 and 9.

Disinfection can be done with 70% disinfectant. Allow to air dry.

Note:

- Never use abrasive and/or harsh detergents or solvents, including thinner, dry cleaning fluid such as thinner, benzene, and the like.
- Do not use a pressure washer to clean the wheelchair.

Please note:

- Tightening the screws with excessive force may damage the wheelchair.

Pneumatic tyres (option)

Check the tire pressure and the condition of the tires (for aging) regularly. In order to drive the wheelchair comfortably, i.e. with the lowest possible rolling and turning resistance, it is important that the tires are properly inflated. When the tire pressure is low, the wheelchair drives and maneuvers more heavily. To inflate the tires, use a valve nipple (e.g. a bicycle pump). The maximum tire pressure is indicated on the tires. Life & Mobility recommends a maximum tire pressure of 6 bar. In case of a flat tire, replace the tire.

Warning:

- Under-inflated tyres have a negative impact on braking power and driving characteristics.

Replacing Quick-Release axles, wheels, or axle sleeves

Note:

Make sure the QR axle is properly adjusted. The wheel must be mounted without any axial play, but at the same time, the safety balls must be securely hooked behind the axle bushing. Check this by pulling the wheels outward; the wheel should not come loose from the frame.

Monthly

Check the functioning and/or damage/signs of wear of the following parts on a monthly basis and replace them if necessary (replace only with original Life & Mobility spare parts):

- Brake
- Tyres (tread and condition of the tyre)
- Seat, backrest, seat mat
- Front wheel, rear wheel and its bearings

Every Quarter

One-Handed Drive (Option)

Regularly lubricate the QR axle with a universal grease designed for high-pressure loaded bearings that operates in wet, dusty, and/or dry conditions, such as Total MULTIS EP 2. Lubricate it in such a way that the axle can rotate smoothly in the axle housing during use of the wheelchair. Assess whether the lubrication interval needs to be shortened due to usage, for example, if the wheels are removed and reattached to the wheelchair more frequently than average.

Major maintenance (annual)

Your dealer will be happy to carry out the major maintenance of the wheelchair for you. This way you can be sure that you are making optimal use of your wheelchair. You can also contact your dealer for the replacement of parts.

Dealer Servicing

During a service, the following points should be checked:

- During a service, the following points should be checked:
- Check parts for breakages. If a break is found, the part must be replaced.
- Check parts for rust. If necessary, touch up with metal paint or replace.
- Check moving parts for smooth movement. If backlash adversely affects operation, modify or replace parts.
- Check bolted connections for tightening torque and, if necessary, apply the overview below; (unless otherwise stated in the manual or instructional videos)
 - M5 screws -> 6 Nm
 - M6 screws -> 10 Nm
 - M8, M10 and M12 screws -> 25 Nm

Reuse

This product is suitable for reuse by another user. The local dealer can prepare the wheelchair for the next user. Cleaning and disinfection must take place before reuse. Parts should be inspected for wear or breakage and replaced if necessary. All necessary adjustments to dimensions and functionality can be made according to the instructions in the manual.

Recalls

Recalls are listed on the website www.life-mobility.com or you will be personally informed about it by your dealer.

Technical Issues

In the unlikely event of immediate technical problems with your wheelchair, we recommend that you contact your dealer. The dealer will assess whether they can repair the wheelchair or whether this needs to be returned to the manufacturer for repair. You may be able to solve some problems yourself.

Adjustable parts such as push handles, armrests, leg rests, etc. are fixed.

- Due to excessive forces on the above parts, the parts may have clamped themselves in the frame. Carefully tap the parts with a plastic mallet.

Allen screws for the settings cannot be loosened

- It may be that the Allen screws have been over-tightened. Try inserting a tube over the long side of the Allen wrench to obtain a longer swivel arm. Now try to loosen the screw.

Please note:

- Do not do this to tighten the screws. This can cause the screw to be over-tightened and cause damage.
- The head of the Allen screw is damaged. Try gripping the screw with a pipe wrench. If this does not work, consult your dealer.

The brakes have insufficient braking power or do not brake with the same braking force on both sides

- You need to readjust the brakes, go to the description on pg. 24.
- Check the tire pressure regularly and inflate the tires if necessary.

Guarantee

Warranty conditions

Life & Mobility's products are carefully assembled and checked before they leave the factory. If it turns out that a product does not meet expectations, you can contact the party from which you purchased the product.

Life & Mobility respects the laws and regulations applicable in the country of purchase to protect consumers.

Defects caused by: non-observance of the instructions for use and maintenance, other than normal use, wear and tear, negligence, overloading, accident by third parties, used, non-original parts and defects caused by the product are not covered by the warranty.

This warranty is void of any other warranty; provided for by law, or communicated orally, except for what is guaranteed in writing by Life & Mobility.

Guarantees are only valid within the EU.

Used wheelchairs/environment

Product disposal and packaging

From an environmental point of view, we recommend that the packaging and the wheelchair be disposed of separately as much as possible so that the materials used can be recycled efficiently.

You can return your used wheelchair to your wheelchair dealer.

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