NEWTON

Newton Curve 1

Newton Curve 2

Newton Curve 3


Newton Curve 3 Deep

Scan the QR code to access all instructional videos on adjusting and using the Newton Curve backrests: Curve 1, Curve 2, Curve 3 and Curve 3 Deep.

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 Life & Mobility

Contents

This manual has been prepared with the utmost care and is based on information known to Life & Mobility at the time of publication.

Information for people with visual impairments can be requested from Life & Mobility. Contact details can be found on the last page.

The information in this manual applies to the standard version of the product. Life & Mobility accepts no liability for damages that may arise from modifications or adjustments to the product after its initial sale. Life & Mobility also cannot be held liable for damages caused by work performed by third parties.

The information in this guide may not be reproduced for commercial purposes and/or published in any form without the prior written permission of Life & Mobility.

According to the Trade Names Act, the trade names and trademarks used by Life & Mobility may not be considered freely available.

No rights can be derived from this text. Subject to change.

© Life & Mobility – All rights reserved

Foreword	4
Product description	7
Technical specifications	10
Safety instructions	12
Settings and customizations	14
Maintenance	20
Guarantee	22
Used wheelchairs/Environment	23

Foreword

These instructions apply to the following Newton Curve backrests: Curve 1, Curve 2, Curve 3, and Curve 3 Deep. Newton Curve backrests are suitable for both indoor and outdoor use. The maximum user weight is 160 kg. This is indicated on the CE sticker, which is located on the pelvic support under the cover.

The Newton backrests Curve 1, Curve 2, Curve 3 and Curve 3 Deep meet the relevant requirements of the wheelchair standard and the MDR (Medical Device Regulation).

Intended purpose:

The backrest of the wheelchair, designed to provide optimal back support and comfort to the wheelchair user.

For the Newton Curve backrests, the following scenarios can reasonably be described that do not correspond to the intended use as described above:

- Failure to heed the warnings in the instructions for use.
- Exceeding the weight and dimension limits specified in the instructions for use.
- Failure to follow the instructions for safe use.

Use

The Newton Curve backrests can be attached to various makes and models of wheelchairs, both electric and manual, using the included mounting hardware. When properly adjusted and positioned, the Newton Curve backrest provides optimal back support for the wheelchair user. Please follow the instructions and step-by-step guide (scan the QR code).

Intended users

The intended users are people with limited mobility, aged between 6 and 100+, and weighing up to 160 kg.

Indications

- Necessary adjustments to a wheelchair in cases of significant or substantial impairment of sitting due to functional and/or structural damage to the trunk or trunk and, if applicable, neck muscles (e.g., due to neurological/neuromuscular diseases, spinal deformities) with postural abnormalities - To maintain the most physiological sitting position possible and ensure stable sitting.

Contraindications

- The Newton Curve backrest must not be used or prescribed by persons who exceed the specified maximum user weight.

Important:

- Read these instructions carefully before using the Newton Curve backrest.
- Failure to comply with safety regulations, instructions, warnings and maintenance recommendations may impair the function of the Newton Curve backrest and may result in injury or damage to the wheelchair or its surroundings.

Tip:

- Sections in this document that are relevant to your safety are marked with a safety symbol.

Safe transport

The Newton Curve backrests have been dynamically tested for use in a motor vehicle. Follow the installation, use, and maintenance instructions in this manual carefully.


Product description

The instructions for your wheelchair product contain information on safe transport in a taxi. Please follow these instructions.

Important:

- If possible, the user should always switch to a fixed seat in the taxi bus.

Universal clamp for back tubes

If you order a Newton backrest from Life & Mobility for mounting on a wheelchair not manufactured by Life & Mobility, this set includes:

0. A Newton Curve backrest of your choice (Curve 1, Curve 2, Curve 3, Curve 3 Deep)
1. Newton Universal back tube mounting clamp (left and right), including Allen key for all adjustments
2. Newton backrest mount
3. Newton backrest insert
4. Newton backrest fitting, suitable for the seat and back width of the wheelchair

The Newton seating system is:

- High-quality back support thanks to 3 adjustable support surfaces.
- Multi-adjustable, without additional parts or tools.
- Easy to adjust while the customer is in the wheelchair, for optimal comfort.
- Compatible with almost all wheelchairs thanks to the universal clamp, which can be easily mounted on all back tubes.

Life & Mobility offers a suitable Curve backrest for every user, regardless of their disability, limitation, or needs. The functional requirements that each backrest fulfills are listed in the table below.

Compatibility

The Newton backrests can be mounted on the following Life & Mobility products: Roxx, Canto Nxt, Neox, and Morgan. They can also be mounted on wheelchairs from other manufacturers, such as the Meyra Netti, Meyra Avanti, Sunrise Rubix, Sunrise Life I, and many more.

Wheelchair dimensions

- The universal Newton Curve backrest mounting kit fits back tubes with a diameter of 19-28 mm and various back tube shapes (not just rectangular back tubes).
- By measuring the distance between the back tubes, you can determine the required Newton back width.

Further information and the latest news on the Newton seating system can be found on the Life & Mobility website at www.life-mobility.com

Universal clamp with QR code

Newton Curve 1 – Front

Newton Curve 1 – Reverse Side

Newton Curve 2 – Front

Newton Curve 2 – Reverse Side

Newton Curve 3 – Front

Newton Curve 3 – Back

Newton Curve 3 Deep – Front

Newton Curve 3 Deep – Back




Technical specifications

Newton Curve Backrests	Curve 1	Curve 2	Curve 3	Curve 3 Deep
Back height (H)	45 cm	45 cm	45 cm	50 cm
Back width (B) ¹	37-57 cm	37-57 cm	37-57 cm	37-57 cm
Contour depth of backrest (T) ²	4 cm	7 cm	15 cm	15 - 25 cm
Wheelchair width	34-65 cm	34-65 cm	34-65 cm	34-65 cm
Angle of inclination	+/- 40°	+/- 40°	+/- 40°	+/- 40°
Weight backrest	From 2.5 kg	From 2.5 kg	From 2.5 kg	From 2.5 kg
Max. user weight	160 kg	160 kg	160 kg	160 kg

Functionality

Adjustable pelvic support	•	•	•	•
Adjustable thoracic support		•	•	•
Lateral support (width adjustment)		•	•	•
Lateral support (height adjustment)			•	•
Adjustable angle of lateral support				•

1. The backrest width can be 2 sizes smaller than the seat width of the wheelchair.
 2. 10 cm deeper with optional side pads

Life & Mobility recommends using the Newton Curve backrest at an ambient temperature between -5 °C and +40 °C.

Service and technical support

For questions regarding specific adjustments, maintenance, or repairs, please contact your dealer. Always have the UDI number ready. You can find it on the CE sticker on the wheelchair, located on the upper frame tube at the front left, under the seat cushion.

Expected lifespan

Life & Mobility specifies a lifespan of 5 years for this product. This means that the product will be repairable, safe, and usable during this period.

The above statement only applies if the user uses the product as intended by Life & Mobility, with CE marking and including all original accessories.

Safety instructions

Life & Mobility accepts no liability for damages or injuries caused by negligence or failure to observe safety instructions.

Depending on the specific circumstances or accessories used, additional safety instructions may be required.

Always inform Life & Mobility and the relevant authority about serious incidents.

General

- Never alter the design of the backrest.
- Use only genuine Life & Mobility parts and accessories. These are tested and safe.
- Do not add any extra weight to the backrest.
- Life & Mobility recommends consulting a doctor or therapist with experience in sitting and positioning to determine if the use of a Newton Curve backrest is appropriate.

Assembly

The Newton Curve backrest may only be installed by the following persons:

- An authorized Life & Mobility dealer.
- A qualified and approved healthcare professional who has completed appropriate training at Life & Mobility or an authorized dealer.
- Before using the wheelchair, check that all parts of the backrest are properly mounted and securely in place.

Warning:

- *Installing a Newton Curve backrest on a wheelchair can change the wheelchair's center of gravity. This creates a risk of the wheelchair tipping backward, which can lead to injuries.*

- *Always check whether precautions need to be taken to prevent the wheelchair from tipping over (e.g., anti-tip supports), or whether adapted fittings need to be fitted for people with an amputation to increase the stability of the wheelchair.*

Use

Warning:

- *Do not use the backrest to push or lift the wheelchair. Improper use can cause the backrest to unexpectedly detach from the wheelchair.*
- *Do not expose the wheelchair to direct sunlight: The padding can become too hot, which can lead to burns or skin sensitivity.*

Settings and customizations

The backrest should ideally be installed without an occupant in the wheelchair.


Once the backrest is installed, minor adjustments can be made while the user is seated in the wheelchair. See the Newton step-by-step plan (scan the QR code).

Required tools


Use a 5 mm Allen key to install and adjust the Newton Curve backrest. This is included in the universal Newton mounting kit.

Mounting the universal mounting adapters

1. Adjust the tilt angle of the back tubes to an angle of approximately 103 degrees to the upper seat frame tube of the wheelchair.
2. Loosen the two Allen screws (A) of the universal mounting adapters enough to allow the clamp to be placed around the back tube.
3. Open the clamp and place it around the backrest tube. The correct individual height position of the backrest will be determined and adjusted later with the person in the wheelchair.
4. Tighten the two Allen screws (A) by hand so that the backrest can later be easily adjusted to the correct individual height.
5. Repeat steps 1 to 3 for the other back tube, mounting the clamps on the left and right at the same height.

Assembly of the Universal mounting adapter

Mounting the backrest fitting on the backrest

Mounting the plug pin on the backrest fitting

Step 1: Adjusting the backrest height

Mounting the backrest fitting on the Newton Curve backrest

1. Attach the backrest fitting to the elongated grey area using the 2 Allen screws (B) in the holes provided for this purpose.
2. Attach the insert (C) to the backrest fitting using the 2 Allen screws (D) and tighten them by hand.
3. Insert the backrest insert (C) into the backrest mounting clamps, whereby the backrest is automatically fixed using the locking tabs (E).

Individual backrest adjustment for the user

These adjustments are made with the person in the wheelchair. These adjustments can be found in the Newton step-by-step plan supplied with the backrest, and the corresponding numbering of the step-by-step plan can also be found on the Newton Curve backrest.

Step 1: Back height (Curve 1, 2, 3 and 3 Deep)

- Loosen both Allen screws (A) of both mounting clamps.
- Adjust the backrest to the correct individual height.
- Position the LT support surface of the backrest (between the two green seams) at the correct height in relation to the user's LT surface (lumbothoracic surface). This begins at the level of the pelvic rim (below the rib cage).
- Tighten both screws (A) of both clamps again to 8 Nm.

Step 2: Backrest angle (Curve 1, 2, 3 and 3 Deep)

- Loosen both screws (D) on the plug pin.
- Adjust the backrest to the correct individual backrest angle.
- Tighten both screws (D) on both sides again.
- Choose a seat surface that matches the backrest angle to prevent slipping.


Step 2: Adjusting the backrest angle

Step 3a: Pelvic support (Curve 1,2,3)

For Curve 1, 2 and 3, follow the steps below; for Curve 3 DEEP, continue at step 3b.

- Remove the protective cap from the Allen screw (F).
- Support the client's individual pelvic position by rotating the pelvic support surface towards the client's pelvis using the Allen screw (F).

Step 3a: Adjusting the pelvic support

Notice:


- Once the support is in place (this can be easily felt between the pelvic support and the patient's buttocks), the support is usually sufficient. Always ensure that no one is pushed forward out of the chair.
- Replace the protective cap on the Allen screw (F).

Step 3b: Pelvic support (Curve 3 Deep)

Position until the pelvis is supported and secure.

Step 3b: Pelvic Support Curve 3 Deep

Step 4a: Adjusting the thoracic support

Step 4a: Thoracic support (Curve 2 and 3)

In order to properly support the individual curvature of the back in the thoracic area, the height (Allen screw G) and the angle (Allen screw H) of the thoracic support are adjusted in this step.

- Set the height (G) to the starting point where the client's kyphotic curvature begins.
- Adjust the angle (H) according to the depth of the kyphosis.

Step 4b: Lumbar-thoracic support

Step 4b: Lumbar-thoracic support (Curve 3 Deep)


- Attach and fasten to the upper back.

-
-
-

Step 5: Lumbar guidance (Curve 2 and 3)

The BOA Fit System ensures that the backrest adapts well to the patient's torso shape at the level of the LT.

- Pull out both BOA buttons (I) to fully relax the enveloping straps.
- Press a BOA button (I) and turn it clockwise to tighten the straps.

Step 5: Adjusting the lumbar support

Danger:

The lumbar support can be adjusted independently on the left and right sides.

Step 6a: Adjusting the lateral support

Step 6a: lateral support

- The height of the lateral support can be adjusted independently on the left and right sides using the two Allen screws (J) on both sides of the backrest.

For additional lateral support, the optional pads can be mounted:

- Loosen the two screws (K) from the fitting of the Newton backrest.
- Place the pads on the backrest and reattach the two screws (K).
- The width position of the pads can be adjusted by slightly loosening the screws (L), moving the pads into the correct position, and tightening the screws again.

Step 6b: lateral support

Step 6b: Lateral support (Curve 3 Deep)

- Adjust the height of the lateral support independently on the left and right sides.
- Adjust the angle of the lateral support using the green lever.

Step 7: Pressure distribution (Curve 3 Deep)

- For optimal pressure distribution, fine-tune the tensioning straps if necessary.


Step 7: Pressure distribution

Remove

After installation and adjustment, the backrest can also be easily removed.

- To remove, push both quick release levers (E) forward to release the lock. *The open lock is now visible (red area).*
- With both levers pushed forward, the backrest can simply be lifted.

Mounting the plug pin on the backrest fitting

Reattach after removal

- To reattach the backrest to the fittings, hold the pins at the same height as the fittings and press the backrest into position.

Danger:

- *When the backrest is inserted correctly, the insert parts lock into place with a "click".*

Warning:

- *If you are unsure whether you heard the click, pull the backrest upwards and check with the lever to ensure it is properly locked in place.*

Assembly instructions Accessories

Headrest

The square headrest tube can be inserted into the integrated headrest recess in the Newton backrests and secured with the star knob. Using the height indicator block attached to the headrest tube, the headrest can be easily reset to the desired individual height after removal.

Relation

The Curve backrests can optionally be equipped with a removable stretch rubber cover, for example in cases of incontinence, heavy soiling, etc.

Maintenance

Maintenance

Danger:

- *Life & Mobility recommends checking all fastening materials for wear and tear every six months, e.g. for loose screws or broken parts.*
- *Loose screws must be tightened according to the assembly instructions. All fasteners must be tightened according to the specifications in the assembly instructions.*
- *Contact your authorized Life & Mobility dealer immediately to replace any broken parts.*
- *Do not use the system if you have found that any parts are loose or broken.*

cleaning

Cleaning the cover

1. Open the zipper on the underside of the cover.
2. Remove the cover from the back panel.
3. Remove the foam and turn the cover inside out.
4. Close the zipper to prevent the Velcro fastener from getting caught on other materials.
5. Machine washable in warm water (60 °C).
6. Allow to drip dry or tumble dry on low heat.
7. Place the foam back into the cover.
8. Close the zipper.

Notice:

- *Do not clean the cover chemically or in industrial washing or drying machines.*
- *Do not clean the cover by steaming.*
- *Do not bleach.*

Cleaning the Newton back shell and mounting parts

- *The back panel can be wiped with soap and water. No further maintenance is required.*

Guarantee

Warranty conditions

Life & Mobility products are manufactured with care and thoroughly inspected before leaving the factory. Should a product not meet your expectations, please contact the retailer from whom you purchased it.

Life & Mobility complies with the laws and regulations in force in the country of purchase to protect consumers.

The warranty excludes defects caused by the following: failure to observe the operating and maintenance instructions, improper use, wear and tear, negligence, overloading, accidents caused by third parties, use of non-original parts and defects whose cause lies outside the product.

This guarantee supersedes all other guarantees, whether statutory or verbal, except for those guaranteed in writing by Life & Mobility.

Guarantees apply exclusively within the EU.

Used wheelchairs/Environment

Recycling/Reuse

This product is suitable for reuse:

- *Have an authorized dealer check the backrest, mounting hardware, foam, and cover for damage or wear.*
- *Clean the entire product.*

Warning:

- *Risk of injury: Always have the product checked for damage or wear by an authorized dealer before passing it on to another user.*
- *Do not use the product (again) if you notice any damage.*

Disposal of product and packaging

For environmental reasons, we recommend disposing of the packaging and the backrest separately as much as possible so that the materials used can be recycled efficiently.

You can return your used backrest to your wheelchair dealer.

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