

# BLAZER

## SPECIFICATIONS BLAZER

|                                    |  |                                    |   |
|------------------------------------|--|------------------------------------|---|
| Seat width                         | 42 – 50,5 cm                               | Action radius                      | 26 - 36 cm                                  |
| Seat depth                         | 42 – 50 cm                                 | Maximum obstacle height            | 5 cm  |
| Seat height (without seat cushion) | 44 cm                                      | Batteries                          | 2 x 50AH (standard)                         |
| Backrest height                    | 46 / 51 cm                                 | Minimum ground clearance           | 7,3 cm                                      |
| Seat angle                         | 8° – 28°                                   | Total weight (including batteries) | 99,5 kg                                     |
| Backrest angle                     | manual: 100° – 110° - electric: 90° - 125° | Maximum user weight                | 136 kg                                      |
| Armrest height                     | 18 – 26 cm                                 | Turning radius                     | 82,5 cm (without user) / 102 cm (with user) |
| Legrest length                     | 37 – 46 cm                                 | Drive                              | RWD   |
| Total height                       | 96,5 cm                                    | Controls                           | VR2   |
| Total length                       | 109 cm                                     | Motor capacity                     | 2 x 350 W                                   |
| Total width                        | 58 – 60 – 63 cm                            | Crash test safe                    | Yes   |
| Power tilt                         | 8° - 28°                                   | Hood colour                        | Black                                       |
| Maximum speed                      | 10 km/h                                    |                                    |   |



### Life & Mobility B.V.

Logistiekweg 7  
 NL- 7007 CJ Doetinchem  
 T +31 (0)314 - 328 000  
 F +31 (0)314 - 328 001

[www.life-mobility.com](http://www.life-mobility.com)  
[info@life-mobility.com](mailto:info@life-mobility.com)



Highly manoeuvrable  
and easy to handle



Inspired by life

02/2021 – 1687AA

Inspired by life



## BLAZER

With an exceptionally narrow 580 mm wheelbase, an easily foldable backrest, a modern design and a very affordable price level, the Blazer is one of the most compact and budget-friendly power wheelchairs on the market.

### Excellent indoor and outdoor manoeuvrability

This rear wheel drive power wheelchair is highly manoeuvrable and easy to navigate indoors as well as outdoors.

The Blazer is crash test approved for transportation in a motor vehicle in accordance with ISO 7176-19. The compact design and the tool-free folding-forward backrest make this wheelchair ideal for both storage and transportation.

### In brief:

- Affordable and modern looking
- Super manoeuvrable and compact for indoor and outdoor use
- Perfect for transportation: Compact (foldable) and crashtest safe
- Easy transfers by flip-back armrests and swing-away legrests
- Three seatwidth options
- 2-wheel suspension and 9-inch front wheel

### Easy transfer

The adjustable, flip-back armrests and swing away footrests make it easy to get on and off the wheelchair or roll closer to a table. The detachable swing-away footrests allows easy transfer, on and off the wheelchair and detaches with no tools required.



**karma**



Contour Adjustable Headrest



Swing Away Controller Bracket



2-in-1 LED Lighting System



Manual Stepless Recline



3 Frame Color Options



Tie-Down Bracket



Anti-Tippers



Tray Table

### STANDARD:

- Flip-back ergonomic armrest
- 2-wheel suspension and 9-inch front wheel
- Adjustable seat depth (42 to 50 cm)
- Adjustable backrest angle  
manual: 100° – 110° - electric: 90° - 125°
- Adjustable armrest height (18 to 26 cm)
- 2-in-1 Led lighting control system

### OPTIONS:

- Power tilt up to 28 degrees
- Elevating legrest with calf support
- Contoured adjustable headrest
- Swing-away controller bracket
- Handle bar
- Push handle
- Tray table



Adjustable backrest angle  
manual: 100° – 110°  
electric: 90°-125°



Foldable backrest for easy storage and Transportation



Adjustable seat depth  
(42 to 50 cm)



Adjustable armrest height  
(18 to 26 cm)

