

BLAZER



Highly manoeuvrable
and easy to handle

Inspired by life

 Life & Mobility

BLAZER

With an exceptionally narrow 580 mm wheelbase, an easily foldable backrest, a modern design and a very affordable price level, the Blazer is one of the most compact and budget-friendly power wheelchairs on the market.

Excellent indoor and outdoor manoeuvrability

This rear wheel drive power wheelchair is highly manoeuvrable and easy to navigate indoors as well as outdoors. The Blazer is crash test approved for transportation in a motor vehicle in accordance with ISO 7176-19. The compact design and the tool-free folding-forward backrest make this wheelchair ideal for both storage and transportation.

In brief:

- Affordable and modern looking
- Super manoeuvrable and compact for indoor and outdoor use
- Perfect for transportation: Compact (foldable) and crashtest safe
- Easy transfers by flip-back armrests and swing-away legrests
- Three seatwidth options
- 2-wheel suspension and 9-inch front wheel

Easy transfer

The adjustable, flip-back armrests and swing away footrests make it easy to get on and off the wheelchair or roll closer to a table. The detachable swing-away footrests allows easy transfer, on and off the wheelchair and detaches with no tools required.



Contour Adjustable Headrest



Swing Away Controller Bracket



2-in-1 LED Lighting System



Manual Stepless Recline



3 Frame Color Options



Tie-Down Bracket



Anti-Tippers



Tray Table

STANDARD:

- Flip-back ergonomic armrest
- 2-wheel suspension and 9-inch front wheel
- Adjustable seat depth (40 - 50 cm)
- Adjustable backrest angle manual: 100° – 114° - electric: 90° - 122°
- Adjustable armrest height (18 to 26 cm)
- 2-in-1 Led lighting control system



Adjustable backrest angle manual: 100° – 114° electric: 90°-122°



Foldable backrest for easy storage and Transportation



Adjustable seat depth (40 to 50 cm)



Adjustable armrest height (18 to 26 cm)

OPTIONS:

- Power tilt up to 28 degrees
- Elevating legrest with calf support
- Contoured adjustable headrest
- Swing-away controller bracket
- Handle bar
- Push handle
- Tray table



SPECIFICATIONS BLAZER

Seat width	41-49 cm	Action radius	25 - 30 km
Seat depth	40 - 50 cm	Maximum obstacle height	5 cm
Seat height (without seat cushion)	44,5 - 61,5 cm	Batteries	2x40Ah (standard)
Backrest height	46 cm	Minimum ground clearance	7,3 cm
Seat angle	8° – 28°	Total weight (including batteries)	102 kg
Backrest angle	manual: 100° – 114° - electric: 90° - 122°	Maximum user weight	136 kg
Armrest height	18 – 26 cm	Turning radius	82,5 cm (without user) / 102 cm (with user)
Legrest length	37 – 46 cm	Drive	RWD
Total height (without headrest))	97,5 cm	Controls	LiNX
Total length (with legrest)	109 – 116 cm	Motor capacity	2 x 350 W
Total width	58 / 60 cm	Crash test safe	Yes
Power tilt	8° - 28°	Hood colour	Black
Maximum speed	10 km/h		

Life & Mobility B.V.

Logistiekweg 7
NL- 7007 CJ Doetinchem
T +31 (0)314 - 328 000
F +31 (0)314 - 328 001

www.life-mobility.com
info@life-mobility.com



Inspired by life