

PRESTOS



Slightly quicker

Inspired by life

 **Life & Mobility**

Slightly quicker

Do you enjoy getting up a bit more speed when you are out and about? The Presto S, the sportier and faster version of the Presto, takes that in its stride. Thanks to its tight turning circle and powerful motor, you can drive nimbly and quickly through the city or over obstacles. The Presto S is equipped with unique daytime running lights that increase your visibility in traffic. Rubber bumpers all round and ample room for your feet mean you will be relaxed when you arrive at your destination. Its sleek design with a matt and gloss black finish will certainly turn heads!

PRESTO^S


17 km/h



44 km

INDEGO

suspension system


SDS
SAFE
DRIVING
SOLUTIONS

Safety

- To increase your visibility, the indicators are positioned high up and the rear lights are designed so they can also be seen from the side.
- The indicators feature daytime running lights that increase your visibility in traffic.
- Directional headlight with LED bulbs for better vision when cornering.
- You can keep your eyes on the road because the control buttons can be identified by touch and are positioned within thumb's reach.
- When cornering, your speed is automatically restricted to a safe level so you remain in full control.
- The indicators switch off automatically after a turn and can be operated safely without having to let go of the steering wheel.

EDS
ERGONOMIC
DRIVING
SOLUTIONS

Comfortable ride

- Outstanding ride comfort thanks to the unique INDEGO suspension system.
- Clearly legible pictograms and information thanks to a large and well laid-out control panel.
- Optimum sitting position thanks to the multi-adjustable comfortable seat: height, depth, backrest angle, armrest height, width and seat cushion angle are fully adjustable.
- Generous leg and foot room thanks to the ergonomic design.

Features

- Adjustable speed up to 17 km/h.
- Maximum operating range of 44 km*.
- Tight turning circle of 160 cm.
- High load capacity of 160 kg.
- Rubber bumpers all round.
- Removable plastic shopping basket.
- Innovative tubeless tyres with low rolling resistance.


* Under ideal circumstances




Adjustable tiller



Rubber impact bumpers



Ergonomic seat



Daytime running lights