

◆ Stand up again, improve your independence  
and be kind to your body!

◀ ERGO STAND ▶



*S-Ergo* seat system

**PATENTED**  
**PATENTED**  
**PATENTED**

😊 **S-ERGO SEAT**

☹️ **TRADITIONAL SEAT**

◀ KP 80 ▶

◆ The **Ergo Stand** is a highly manoeuvrable standing power wheelchair with *Karma's* unique, patented *S-Ergo* seating system. The front wheel drive *Ergo Stand* is only 62 cm wide and has excellent indoor manoeuvrability performance. Agile control and pleasant driving characteristics during outdoor use is guaranteed. The *Ergo Stand* integrates many adjustable functions to meet individual user requirement. The power recline function provides a wide angle adjustment of 80° – 120° with excellent back support.

- ◆ Power stand
- ◆ Power backrest recline



# Features



• FRONT WHEEL DRIVE



• POWER SIT TO STAND



• POWER RECLINE



• FLIP BACK AND HEIGHT ADJUSTABLE ARMREST



• LIGHTS



• KNEE SUPPORT



• SWING AWAY JOYSTICK



Reduced sliding



Pressure redistribution

• S-ERGO SEAT SYSTEM

## < The medical benefits of standing >

- ◆ Activates circulation system
- ◆ Improves cardio-pulmonary function
- ◆ Helps prevent decubitus sores
- ◆ Improves bowel regularity and function
- ◆ Helps bladder and intestinal motility
- ◆ Reduces contractures and ankylosis

## SPECIFICATION

STANDARD	Speed:	6 – 9 km/h
	Width:	62 cm
	Length:	100 cm
	Seat width:	40 / 44 / 48 cm
	Seat depth:	45 / 48 / 51 cm
	Motors:	2x 300 W
	Batteries:	2x 50 Ah
	Range:	25 km/h
	Max. user weight:	135 Kg
	Wheelchair weight:	83 Kg
	Electronics:	VR2

S-Ergo seat system
Power stand
Power recline backrest
Flip back and height adjustable armrests
Knee support adjustable in height, depth and width
Light kit
Adjustable footplate height with 7 levels
Front wheel drive
Color: Diamond blue
Driving in standing position
± 1cm ± 0,5 kg

OPTIONS	Parallel swing-away joystick	Seat belt
	Head support	Mud guard
	H-shaped seat belt	

Karma Europe reserves the right to modify information herein without further notice.

