

◆ *The most economic power wheelchair, without sacrificing safety and durability.*

## < ERGO TRAVELLER >



◆ The *Ergo Traveller* is Karma's most economic power wheelchair. Its lightweight and compact design guarantees outstanding maneuverability during indoor use. The powerful base also ensures great outdoor performance. The concept behind the design of the *Ergo Traveller* was to make it as light as possible and convenient to transport. The patented *S-Ergo* system promotes pressure redistribution, reduces downward sliding and helps maintain good posture. The *Ergo Traveller* can be taken apart quickly and easy into three easy to handle components. The heaviest part weighs just 24 kg. The quick release battery box has a convenient charging point built into it, which allows the batteries to be charged away from the power chair.

*S-Ergo* seat system


## < KP 10 >



• BATTERY BOX



• EASY DISASSEMBLED INTO 3 PARTS



• FLIP BACK ARMREST

### SPECIFICATION

STANDARD	Speed:	6 km/h
	Width:	62 cm
	Length:	103 cm
	Seat size:	43 x 45 cm
	Motors:	2x 200 W
	Batteries:	2x 22 Ah
	Range:	16 – 20 km/h
	Max. user weight:	115 Kg
	Wheelchair weight without batteries:	25 Kg
	Electronics:	VR2

OPTIONS	Elevating and swing-away footrest
	Attendant joystick kit
	Paralell swing-away joystick

- Ergonomic seat system provides comfortable seating
- Tool-free disassembled into 3 parts for traveling and transportation
- Quick release battery pack
- Washable upholstery
- Height adjustable and flip back armrest
- 2 staged adjustable back angle
- 3 staged height adjustable footrest
- Right or left hand joystick control positions
- Lightweight design, the heaviest part after taking apart is only 24 kg
- ± 1cm ± 0,5 kg

**i** Karma Europe reserves the right to modify information herein without further notice.

