



# Product overview

**ROLLIT**



## Rollit - Basic model

|                  |                       |
|------------------|-----------------------|
| Frame colour     | RAL 1028              |
| Seat width       | 47 cm                 |
| Seat height      | 51 cm                 |
| Type of legrests | Footrest + heel strap |
| Seat             | Standard seat         |
| Backrest         | Standard backrest     |
| Sides            | Folding sides         |



## Specifications

| Abbr. | Measurement points     | Specifications |
|-------|------------------------|----------------|
| SW    | Seat width             | 47 / 52 cm     |
| SD    | Seat depth             | 46 cm          |
| SH    | Seat height            | 51 / 56 cm     |
| BH    | Backrest height        | 47 cm          |
| AH    | Armrest height         | 24 cm          |
| WW    | Width of wheelchair    | 72 cm          |
| HP    | Height of push-handles | 92 / 97 cm     |
| LL    | Lower leg length       | 45 - 55 cm     |
|       | Total weight           | 39 kg          |
|       | Maximum user weight    | 160 kg         |





## Product overview

### Contact details

|               |       |
|---------------|-------|
| Dealer's name | _____ |
| Contact       | _____ |
| Date          | _____ |
| Order number  | _____ |

### Model

|                                                                                                                                                                                              |         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|
| <input checked="" type="checkbox"/> Rollit                                                                                                                                                   | € 1.433 |
| <i>Rollit frame: Rigid internal frame, folding sides with long armrests, 7" semi-solid front tyres, 24" semi-solid rear tyres, upholstery colour Midnight blue and extended wheelbrakes.</i> |         |
| Quantity .....                                                                                                                                                                               |         |

Please tick what is required! The options listed may be subject to an additional charge.

### Frame colour

|                                                    |   |
|----------------------------------------------------|---|
| <input type="checkbox"/> RAL colour - No.1,2 ..... | - |
|----------------------------------------------------|---|

### Seat width

|                                |       |
|--------------------------------|-------|
| <input type="checkbox"/> 47 cm | -     |
| <input type="checkbox"/> 52 cm | € 210 |

### Seat height

|                                |       |
|--------------------------------|-------|
| <input type="checkbox"/> 51 cm | -     |
| <input type="checkbox"/> 56 cm | € 102 |

### Legrest

|                                                                                                           |       |
|-----------------------------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> L <input type="checkbox"/> R Footrest + heel strap                               | € 123 |
| <input type="checkbox"/> L <input type="checkbox"/> R Comfort legrest                                     | € 293 |
| <input type="checkbox"/> L <input type="checkbox"/> R Comfort legrest + heel strap                        | € 304 |
| <input type="checkbox"/> Lower leg length adjustable by means of allen wrench                             | -     |
| <input type="checkbox"/> Lower leg length adjustable by means of wingnut                                  | -     |
| <input type="checkbox"/> L <input type="checkbox"/> R Knee cushion                                        | € 77  |
| <input type="checkbox"/> L <input type="checkbox"/> R Amputation support <sup>3</sup>                     | € 257 |
| <input type="checkbox"/> Length 25 cm                                                                     | -     |
| <input type="checkbox"/> Length 40 cm                                                                     | -     |
| <input type="checkbox"/> Terry cover amputation support ..... pcs                                         | € 49  |
| <input type="checkbox"/> L <input type="checkbox"/> R Plaster cast support <sup>3,4</sup>                 | € 330 |
| <input type="checkbox"/> L <input type="checkbox"/> R Pes equinus cushion as part of plaster cast support | € 132 |
| <input type="checkbox"/> Terry cover pes equinus cushion ..... pcs                                        | € 56  |

### Seat

|                                                          |       |
|----------------------------------------------------------|-------|
| <input type="checkbox"/> Standard seat                   | € 168 |
| <input type="checkbox"/> Plank seat                      | € 179 |
| <input type="checkbox"/> Tuber opening for Plank seat    | € 41  |
| <input type="checkbox"/> Front rim h=5 cm for Plank seat | € 56  |
| <input type="checkbox"/> Front rim h=8m for Plank seat   | € 71  |
| <input type="checkbox"/> Contour+ seat                   | € 459 |

### Backrest

|                                              |       |
|----------------------------------------------|-------|
| <input type="checkbox"/> Standard backrest   | € 173 |
| <input type="checkbox"/> Anatomical backrest | € 184 |

### Sides

|                                                                        |       |
|------------------------------------------------------------------------|-------|
| <input type="checkbox"/> Folding sides                                 | -     |
| <input type="checkbox"/> Folding sides with height-adjustable armrests | € 174 |
| <input type="checkbox"/> Fixing clip set for sides                     | € 21  |

### Headrest

|                                                      |       |
|------------------------------------------------------|-------|
| <input type="checkbox"/> Standard headrest small     | € 199 |
| <input type="checkbox"/> Standard headrest           | € 221 |
| <input checked="" type="checkbox"/> Headrest adapter | € 51  |

### Worktop

|                                                                               |       |
|-------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> Perspex worktop                                      | € 361 |
| <input type="checkbox"/> Trespa worktop                                       | € 235 |
| <input type="checkbox"/> L <input type="checkbox"/> R College Perspex worktop | € 331 |

### Accessories

|                                                                                        |                  |
|----------------------------------------------------------------------------------------|------------------|
| <input type="checkbox"/> Stick-holder                                                  | € 60             |
| <input type="checkbox"/> Anti-tipping support <sup>3</sup>                             | € 198            |
| <input type="checkbox"/> Catheter hook                                                 | € 35             |
| <input type="checkbox"/> Oxygen cylinder holder 12,5 cm <sup>3</sup>                   | € 117            |
| <input type="checkbox"/> Spoke protectors                                              | € 246            |
| <input type="checkbox"/> Spoke protectors with printed logo <sup>5</sup>               | Price on request |
| <input type="checkbox"/> Logo plate on backrest with printed logo <sup>5</sup>         | Price on request |
| <input type="checkbox"/> L <input type="checkbox"/> R Adapter for IV pole <sup>3</sup> | € 103            |
| <input type="checkbox"/> IV pole, adjustable in height                                 | € 196            |
| <input type="checkbox"/> Anti-theft pole <sup>3</sup>                                  | € 109            |
| <input type="checkbox"/> Fitted with fixing tube for logo purposes <sup>5</sup>        | € 25             |
| <input type="checkbox"/> Fitted with 2 infusion hooks                                  | € 41             |
| <input type="checkbox"/> Flag with printed logo <sup>5</sup>                           | Price on request |

### Explanation

1 Next RAL colours are not available: RAL 1026 / 2005 / 2007 / 3024 / 3026 / 6038

2 More than 175 RAL colours available without surcharge. Please note: a small number of colours is subject to a surcharge. Just contact your sales advisor!

3 Only available in black (metal parts)

4 Left + right plaster cast support only available from seat width 47 cm

5 Please note: for a printed logo, we need an .eps file of the logo concerned

For changes to the options specified above, please refer to the supplementary order form for alterations.